AlphaDie MF

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Directions for use

Range of Applications

AlphaDie MF is excellent for fabricating precision models for all types of crown and bridgework as well as for dimensionally stable reproduction of shoulders and partial crown preparations.

AlphaDie Separator must be applied when casting polyether impression, resp., duplicating materials.

Do not use with alginates and hydrocolloid impression materials.

Properties

AlphaDie MF is flowable, accurate in detail and dimensionally stable. It has high edge strength and abrasion resistance and is easy to trim and saw.

Comprises

AlphaDie MF is supplied in glass bottles (200 ml base material, 100 ml hardener) and in a plastic container (400 g filler). A dispensing scoop for the filler, mixing/dispersing vessel and spatula are included in the set.

Mixing AlphaDie MF

Briefly shake the bottles containing both the base material and hardener prior to mixing. Note the coloured marking on the cap. Close the glass bottles carefully immediately after use. **Do not allow AlphaDie MF to come into contact with water** (foam).

- Fill the required amount of base material in the dispensing and mixing container supplied. Then add the correct quantity of filler and thoroughly spatulate the mixture. Add the correct quantity of hardener and spatulate the mixture again thoroughly.

Mixing ratio:

Base: hardener = 2 : 1.

Approx. 15 g AlphaDie MF and filler is required for a full dental arch.

10 g base

1 level measuring scoop of filler (approx. 14 g)

5 g hardener

- After mixing for approx. 30 sec., pour the AlphaDie MF in a thin stream into the cleaned and dried impression. The material remains flowable for approx. 2 min. at 20°C.

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- In order to prevent incorporation of air voids, AlphaDie MF may be hardened for 15 mins., after pouring, in a dry pressure vessel at 2-4 bars.

AlphaDie MF is sufficiently hard after 1 hour that the arch can be trimmed, drilled and ground.

Hazard statements

AlphaDie MF Basis

Contains benzene, C10-13-alkyl derivatives

Danger: Causes skin irritation. Causes serious eye irritation. May be fatal if swallowed and enters airways.

AlphaDie MF Härter

Contains diphenylmethane diisocyanate, isomeres and homologues

Danger: May cause allergy or asthma symptoms or breathing difficulties if inhaled. May cause an allergic skin reaction. Suspected of causing cancer.

Precautionary statements

AlphaDie MF Basis:

IF SWALLOWED: Immediately call a POISON CENTER/doctor. Do NOT induce vomiting. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Take off contaminated clothing and wash it before reuse.

AlphaDie MF Härter:

Wear protective gloves/protective clothing/eye protection/face protection/hearing protection. In case of inadequate ventilation wear respiratory protection.

Additional information on

AlphaDie MF Härter:

Effective 24 August 2023 appropriate training is mandatory prior to industrial or professional use. Contains isocyanates. May produce an allergic reaction.

CAUTION!

Do not overheat AlphaDie MF. Steam-clean AlphaDie MF (max. 15 secs.).

Disposal

Dispose of in a special refuse tip, in accordance with local regulatory requirements.

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