FuturAcryl 2000

Schütz Dental GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany
Telefon: +49 (0) 6003 814-0 · Telefax: +49 (0) 6003 814-900
info@schuetz-dental.de · www.schuetz-dental.de

CE 0297

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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Instructions for use

General product description:

Heat-curing acrylic for the production of acrylic dentures in the injection technique or in the flask technique.

Composition:

- Powder: bead polymer made of polymethyl methacrylate (PMMA), pigments, initiators
- Liquid: methyl methacrylate, tetramethylene dimethacrylate, initiators, stabilizers
- FuturAcryl 2000 is free of cadmium.

Indications:

- For the production of partial and full dentures in the injection technique or in the flask technique.
- For the production of splints
 - occlusal splints
 - reflex / anti-grinding splints

Contraindications / interactions with other agents:

Do not use in case of a known allergy against one of the components.

Application / processing:

Preparatory steps:

Produce plaster model and wax-up as usual. The prefabricated teeth set up in the wax are fixed with plaster or silicone (stent / flask). Subsequently, the wax is eliminated as usual.

Roughen the acrylic teeth, provide them with undercuts and coat them with monomer or HS Cross-Liquid (Schütz Dental) to ensure that they are securely retained. The surfaces of the plaster models are insulated with the alginate separating agent FuturaSep Plus (Schütz Dental) according to the manufacturer's instructions.

Mixing ratio:

2.5 parts powder by weight, 1 part liquid by weight.

If a longer working range is required, we recommend a mixing ratio of 2 parts powder by weight and 1 part liquid by weight.

Veins can be added, if desired. Please note the separate manufacturer's information and instructions.

Production of acrylic dentures:

After adding the powder to the liquid, stir FuturAcryl 2000 thoroughly with a spatula for 30 seconds in a mixing cup to create a homogenous mixture.

After allowing a doughing time of ca. 10 min., the material is ready for use.

The working range (21 °C / 69 °F) is: ca. 1 hour (with a mixing ratio of 2.5 : 1) or 2 hours (with a mixing ratio of 2 : 1). After the doughing time, knead the acrylic dough for ca. 30 seconds.

Once FuturAcryl 2000 becomes viscoelastic, pack it into the prepared flask and press. Hold the pressure for ca. 15 minutes. Secure the flask in a clamp for polymerisation.

Polymerize FuturAcryl 2000 in the closed flask with the clamp for 20 minutes in boiling water (or in a pressure pot at 95 °C / 203 °F). Allow the flask to cool to room temperature in the clamp and then deflask in the usual way.

Processing in the UNIPRESS system:

Prepare flasks and plaster surfaces as described in the system's instructions for use. Mix FuturAcryl 2000 as specified in the pressing technique. Pack the material into the cylinder and inject. The polymerization procedure corresponds to the pressing technique.

Finishing the material:

Finish the material with suitable tools. Pre-polish with sand paper or pumice.

Polishing:

Polish the material with goat hair brushes, polishing paste and soft woolen wheels. A careful surface treatment and polish are a precondition for an optimal result.

Repairs:

Cold curing resins (e.g. FuturaGen, FuturaPress LT, FuturaSelf or FuturaPress N, Schütz Dental) can be used for repairs. Roughen the fracture surfaces, clean them and coat them with monomer before applying the acrylic dough. The polymerization time is 15 min. in a pressure polymerization

system, at a pressure of 2 - 5 bar and a temperature of 45 °C (113 °F). Please observe the manufacturer's separate instructions for use.

Production of splints:

Preparatory steps:

Impression and production of a plaster model.

Blocking undercuts with suitable materials. Carve splint in wax as desired.

Use plaster to produce a counter model and remove wax.

Insulate plaster model with suitable materials (e.g. FuturaSep, Schütz Dental) according to the manufacturer's instructions.

Mixing ratio:

2.5 parts powder by weight, 1 part liquid by weight.

If a longer working range is required, we recommend a mixing ratio of 2 parts powder by weight and 1 part liquid by weight.

Production of splints:

After adding the powder to the liquid, stir FuturAcryl 2000 thoroughly with a spatula for 30 seconds in a mixing cup to create a homogenous mixture.

After allowing a doughing time of ca. 10 min., the material is ready for use.

Apply the material to the model and close the flask. Secure the flask in a clamp and press. Hold the pressure for ca. 15 min.

The flask is secured in a clamp for polymerization. Polymerize FuturAcryl 2000 in the closed flask with the clamp for 20 minutes in boiling water (or in a pressure pot at 95 °C / 203 °F). Allow the flask to cool to room temperature in the clamp and then deflask in the usual way.

Finishing the material:

After polymerization, the splint is lifted from the model and it is finished with suitable tools.

Silicone polishers (wheel, lentil, roll) as well as hard metal cutters or diamond grinding tools are suited for finishing the material.

Prepolish with sandpaper or pumice.

Polishing:

Polish the material with goat hair brushes, polishing paste and soft woolen wheels. A careful surface treatment and polish are a precondition for an optimal result.

Repairs:

Cold curing resins (e.g. FuturaGen, FuturaPress LT, FuturaSelf or FuturaPress N, Schütz Dental) can be used for repairs. Roughen the fracture surfaces, clean them and coat them with monomer before applying the acrylic dough. The polymerization time is 15 min. in a pressure polymerization system, at a pressure of 2 - 5 bar and a temperature of 45 °C (113 °F).

Please observe the manufacturer's separate instructions for use.

Troubleshooting / FAQ Liste:

Problem	Cause	Remedy
acrylic does not cure properly	incorrect mixing ratio	<ul style="list-style-type: none"> – Gobserve instructions for use – observe mixing ratio
acrylic shows a whitish tint along the bottom	<ul style="list-style-type: none"> – incorrect or improper insulation – model was not soaked in water 	<ul style="list-style-type: none"> – insulate properly – check insulation agent, if necessary – soak model in lukewarm water for ca. 10 min.
acrylic forms bubbles along the bottom	– model was not soaked in water	soak model in lukewarm water for ca. 10 min.
acrylic surface white / spongy	too much time passed until the acrylic was injected or until the acrylic was put in the pressure pot	observe processing times
color deviations of repairs	incorrect mixing ratio	observe mixing ratio
acrylic fractures when lifted	improperly insulated	check insulation
acrylic / prosthesis shows whitish streaks	<ul style="list-style-type: none"> – improperly mixed – powder content too high 	<ul style="list-style-type: none"> – mix acrylic dough thoroughly with a spatula – observe mixing ratio
splint fractures	design too thin / unsuitable	ensure a sufficient layer thickness

Guarantee:

Our processing instructions, regardless of whether they are issued verbally, in writing or in the form of a practical demonstration, are based upon our own experience and can therefore be regarded solely as guidelines. Our products are subject to continuous development. We therefore reserve the right to make possible changes.

Storage:

Storage temperature at 10 - 25 °C / 50 - 77 °F.
Close container carefully.

Shelf-life:

The maximum shelf-life is printed on the label of each packaging. Do not use after the expiration date.

Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn

about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-indications / interactions:

If a patient has known hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth when using the product.

Note:

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Disposal:

Powder can be disposed of as domestic waste. Liquid must be disposed of as hazardous waste.

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