# FuturaGen

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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## Instructions for use

### Composition:

- Powder: pearl polymer made of polymethyl methacrylate (PMMA), pigments, initiators
- Liquid: methyl methacrylate (MMA), tetra-methylene dimethacrylate, initiators, stabilizers
- FuturaGen is free of cadmium

### Indications:

- For the production of partial and full dentures in the injection technique, in the flow technique and in the flask technique.
- For the production of the acrylic part of implant-supported restorations and on frameworks for dentures.
- For repairs and relinings of existing restorations.
- For the production of splints
  - occlusal splints
  - reflex / anti-grinding splints

### Contraindications / interactions with other agents:

Do not use in case of a known allergy against one of the components.

### Advantages of the combination of FuturaGen and the UNIPRESS system:

- optimal reproduction of the occlusion
- excellent fit of partial and full dentures
- flowable – thus no contact with the hands
- economic work

### Application / processing:

#### Preparatory steps:

Produce plaster model and wax-up as usual. The prefabricated teeth set up in the wax are fixed with plaster or silicone (stent / flask). Subsequently, the wax is eliminated. Roughen the acrylic teeth, provide them with undercuts and coat them with monomer or HS Cross-Liquid (Schütz Dental) to ensure that they are securely retained. The surfaces of the plaster models are insulated with the alginate separating agent FuturaSep Plus (Schütz Dental) according to the manufacturer's instructions.

#### Mixing ratio:

For the UNIPRESS technique, we recommend a mixing ratio of 14 weight parts powder, 6 weight parts liquid.

56 ml powder : 16 ml liquid (volume) or 35 g powder : 15 g liquid (weight) are required for an average partial denture.

Veins can be added, if desired. Please note the separate manufacturer's information and instructions.

### **Production of partial and full dentures as well as of implant-supported restorations and framework completions in the injection system UNIPRESS (please observe the separate instructions) (in case of other injection systems, please observe the manufacturer's instructions):**

Mix powder and liquid (according to the mixing ratio). After allowing a doughing time of 15 sek., stir FuturaGen with a spatula in a porcelain mixing bowl to produce a homogeneous mixture and fill into the UNIPRESS cylinder.

As soon as the material has reached a viscoelastic consistency, inject it into the prepared UNIPRESS flask.

These work steps must be completed after 8 min. at the very latest. Polymerization takes place under a defined pressure build-up inside the plate-and-spring system inside of UNIPRESS. The FuturaGen dentures can be deflasked after 20 - 30 min.

### **Production of partial and full dentures as well as of implant-supported restorations and framework completions in the flow technique or in the flask technique:**

Mix powder and liquid (according to the mixing ratio) and fill the mold as usual.

At room temperature (21 °C / 69 °F), FuturaGen is pourable for ca. 3 min. After 5 - 6 min., the material enters the plastic phase and after 7 - 8 min., at the very latest, the restoration must be placed in the pressure pot.

When using **FuturaGen LF Liquid**, the acrylic is pourable for ca. 3 - 4 min., it will enter the plastic phase after 7 - 8 min. and the restoration must be placed in the pressure pot after 9 - 10 min., at the very latest.

Polymerize FuturaGen in a pressure polymerization system for 15 min. at a pressure of 2 - 5 bar and a water temperature of 45 °C / 113 °F.

### **Repairs and relining of existing appliances:**

Produce a mold or a stent, depending on the construction.

Clean the acrylic surfaces to be repaired, roughen them and smooth out the edges. Wet with

monomer or with HS Cross-Liquid (Schütz Dental). Mix acrylic according to the mixing ratio and let the acrylic flow in.

Polymerize FuturaGen in a pressure polymerization system (e. g. Futuramat by Schütz Dental) for 15 min. at a pressure of 2 - 5 bar and a water temperature of 45 °C / 113 °F.

**Finishing the material:**

After polymerization, the restoration is lifted from the model and it is finished with suitable tools.

Silicone polishers (wheel, lentil, roll) as well as hard metal cutters or diamond grinding tools are suited for finishing the material.

Prepolish with sandpaper or pumice.

**Polishing:**

Polish the material with goat hair brushes, polishing paste and soft woolen wheels. A careful surface treatment and polish are a precondition for an optimal result.

**Production of splints:**

**Preparatory steps:**

Production of a plaster model.

Blocking undercuts with suitable materials.

Carve splint in wax as desired.

Produce impression, fit it with sprue canals and remove wax.

Insulate plaster model with suitable materials (e. g. FuturaSep, Schütz Dental).

**Processing:**

Mix powder and liquid (according to the mixing ratio) and fill the mold as usual.

These work steps must be completed after 7 - 8 min. at the very latest.

For a working range extended by 1 - 2 min., use the FuturaGen LF Liquid.

Polymerize FuturaGen in a pressure polymerization system (e. g. Futuramat, Schütz Dental) for 15 min. at a pressure of 2 - 5 bar and a water temperature of 45 °C / 113 °F.

**Finishing the material:**

After polymerization, the splint is lifted from the model and it is finished with suitable tools.

Silicone polishers (wheel, lentil, roll) as well as hard metal cutters or diamond grinding tools are suited for finishing the material.

Prepolish with sandpaper or pumice.

**Polishing:**

Polish the material with goat hair brushes, polishing paste and soft woolen wheels. A careful surface treatment and polish are a precondition for an optimal result.

**Repairs:**

Produce a die or a stent, depending on the construction.

Clean the acrylic surfaces to be repaired, roughen them and smooth out the edges. Wet with monomer or with HS Cross-Liquid (Schütz Dental). Mix acrylic according to the mixing ratio and let the acrylic flow in.

Polymerize FuturaGen in a pressure polymerization system (e. g. Futuramat, Schütz Dental) for 15 min. at a pressure of 2 - 5 bar and a water temperature of 45 °C / 113 °F.

## Troubleshooting / FAQ Liste

<b>Problem</b>	<b>Cause</b>	<b>Remed</b>
acrylic does not cure properly	incorrect mixing ratio	<ul style="list-style-type: none"> <li>- observe instructions for use</li> <li>- observe mixing ratio</li> </ul>
acrylic shows a whitish tint along the bottom	<ul style="list-style-type: none"> <li>- incorrect or improper insulation</li> <li>- model was not soaked in water</li> </ul>	<ul style="list-style-type: none"> <li>- insulate properly</li> <li>- check insulation agent, if necessary</li> <li>- soak model in lukewarm water for ca. 10 min.</li> </ul>
acrylic forms bubbles along the bottom	- model was not soaked in water	soak model in lukewarm water for ca. 10 min.
acrylic surface white / spongy	too much time passed until the acrylic was injected or until the acrylic was put in the pressure pot	observe processing times
color deviations of repairs	incorrect mixing ratio	observe mixing ratio
acrylic fractures when lifted	improperly insulated	check insulation
acrylic / prosthesis shows whitish streaks	<ul style="list-style-type: none"> <li>- improperly mixed</li> <li>- powder content too high</li> </ul>	<ul style="list-style-type: none"> <li>- mix acrylic dough thoroughly with a spatula</li> <li>- observe mixing ratio</li> </ul>
acrylic polymerizes too quickly	<ul style="list-style-type: none"> <li>- temperatures too high</li> <li>- work steps too long</li> </ul>	<ul style="list-style-type: none"> <li>- observe temperatures</li> <li>- observe working range</li> </ul>
splint fractures	design too thin / unsuitable	ensure a sufficient layer thickness (min. 2 - 3 mm)

**Guarantee:**

Our processing instructions, regardless of whether they are issued verbally, in writing or in the form of a practical demonstration, are based upon our own experience and can therefore be regarded solely as guidelines. Our products are subject to continuous development. We therefore reserve the right to make possible changes.

**General notes:**

For deviations from 21 °C room temperature: Higher temperatures reduce the above times, lower temperatures increase them. To further reduce the risk of intolerance reactions, we recommend storing the orthodontic appliance in lukewarm water for at least 12 hours before insertion.

**Storage:**

Storage temperature at 10 - 25 °C / 50 - 77 °F. Close container carefully.

**Shelf-life:**

The maximum shelf-life is printed on the label of each packaging. Do not use after the expiration date.

**Side-effects:**

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

**Contra-indications / interactions:**

If a patient has known hyper-sensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, the composition of the medical device we supply is available upon request. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth before using the product.

**Note:**

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

**Disposal:**

Powder can be disposed of as domestic waste. Liquid must be disposed of as hazardous waste.

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