# Micronium Exclusiv

# Micronium N10

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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### Instruction for use

#### Intended use

Micronium Exklusiv / Micronium N10 is a spring-hard non-precious metal CoCr alloy type 5 for the fabrication of removable dentures. This product is beryllium- and nickel-free, which means that it contains less than or equal to 0.1% beryllium and nickel by mass. **Processing by professional users only!**

#### General guidelines for handling

This instruction for use includes important processing steps and recommendations for Micronium Exklusiv / Micronium N10.

#### Indication

For the production of model cast works for removable restorations.

#### Design

The design has to be made on an investment model under consideration of the standard rules of designs for partial dentures. The modelation should be done with wax that fire without leaving residues.

#### Sprue design

The sprues should have a diameter of 3.5 mm with a length of 20 - 30 mm. Massive areas should additional have a melt reservoir (Ø 5 mm). The cast hopper should be 4 mm above the highest point of the partial denture wax design.

#### Melting and casting

Micronium Exklusiv / Micronium N10 should be melted in a **ceramic crucible**. **Please do not use graphite crucibles and no flux!** Avoid the overheating of the melt. Prevent multiple casts of melt bottoms. The chemical and mechanical properties can only be guaranteed for new material. Melting with open flame (acetylene/ oxygen) and inductive melting: Once the cylinders are melted and a so-called shadow falls across the molten metal, before the oxide skin begins to split, start the casting. **Maximum temperature for casting: 1490 °C.**

#### Devesting

Let the muffle cool down to room temperature (ca. 20 °C), do not quench with water. Put the cooled muffle into water to avoid dust generation during the devesting. Sandblast the surface with 250 µm of aluminium oxide with 3 - 4 bar. Clean

the partial denture base with a steam cleaner.

#### Elaboration

The framework can be elaborated with standard carbide cutters or aluminium oxide stones and rubbers, look for smooth transitions.

#### Soldering / Laser welding

Micronium Exklusiv / Micronium N10 can be soldered with all suitable solder. Micro-nium Exklusiv / Micronium N10 parts should not be soldered with gold or palladium solder. Micronium Exklusiv / Micronium N10 3is also ideally suitable for laser-welding.

#### Handling conditions / Safety

**Metal dust is harmful to health. Use when grinding and sandblasting dust extraction and respirator with filter FFP3 – DIN EN 149.**

#### Contraindications and side effects

If the instructions are observed during the production processes, incompatibilities with non-precious dental alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the contra-indications and side effects. Any serious incident that involve the product must be reported to the manufacturer and the competent authority in the accorded country.

#### Single-use

The chemical and mechanical properties can only be guaranteed for new material.

#### Disposal Instructions

Consult the material safety data sheets or national regulations for disposal. Dispose of Micronium Exklusiv / Micronium N10 residues and dust in an environmentally friendly manner. Grinding dust must not enter groundwater, water bodies or sewers. Address waste exchanges for recycling.

#### Storage conditions

Temperature, humidity or light has no effect on the product properties.

#### Quantity

Please consider the label on the package.

*Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.*

#### Micronium Exklusiv

#### Chemical Composition

Co <span> </span> %	Cr <span> </span> %	Mo <span> </span> %	Fe <span> </span> %
<b>64,2</b>	<b>28,5</b>	<b>5,2</b>	<b>0,5</b>
Si <span> </span> %	Mn <span> </span> %	C <span> </span> %	Others <span> </span> %
<b>0,9</b>	<b>0,3</b>	<b>0,3</b>	<b>&lt; 0,1</b>

Determination, extent and validity of the chemical composition according to DIN EN 10 204 – 3.1

#### Typical technical Data

Yield strength 0,2%	<b>540 MPa</b>
Elongation	<b>5%</b>
Tensile strength	<b>770 MPa</b>
E-module	<b>230 MPa</b>
Density	<b>8,3 g / cm³</b>
Corrosion resistance	<b>&lt; 200 µg / cm²</b>
Tarnish resistance	<b>yes</b>
Hartness	<b>390 HV 10/30</b>
Melting range (Solidus / Liquidus)	<b>1340 °C / 1410 °C</b>
Casting temperature	<b>1490 °C</b>

#### Micronium N10

#### Chemical Composition

Co <span> </span> %	Cr <span> </span> %	Mo <span> </span> %	Si <span> </span> %
<b>67,4</b>	<b>26</b>	<b>5</b>	<b>0,5</b>
Fe <span> </span> %	C <span> </span> %	Mn <span> </span> %	Others <span> </span> %
<b>0,45</b>	<b>0,4</b>	<b>0,3</b>	<b>&lt; 0,1</b>

Determination, extent and validity of the chemical composition according to DIN EN 10 204 – 3.1

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