# Nicor

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product. As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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Mandler 05/2024 - 500

**Instruction for use****Intended use**

Nicor is a beryllium free nickel based bonding alloy for the production of crowns and bridges. Nicor should not be used in patients with a known sensitization to nickel.

**General guidelines for handling**

This instruction for use includes important processing steps and recommendations for Nicor.

**Indication**

For the production of crowns and bridges for the ceramic veneering.

**Modelation**

The modelation should be done with wax that fire without leaving residues under consideration of the standard rules of designs for dental technicians. The frame has to be designed in an anatomical reduced form. The wall thickness should be at a minimum of 0.4 mm to secure the flow out of the melt. Consider a sufficient connector (6-9 mm<sup>2</sup>). Avoid sharp edges and undercuts.

**Sprue design**

We recommend the design of the sprue with a bar. The horizontal sprue should have Ø 4-5 mm, the sprue to the restoration should have Ø 3 mm. Single crowns should be directly connected with a sprue of Ø 4 mm with a length of 15-20 mm. Connect the sprue on massive areas e. g. palatal and avoid the center of the muffle.

**Melting and casting**

Nicor should be melted in a ceramic crucible. Please do not use graphite crucibles and no flux! Avoid the overheating of the melt. Prevent multiple casts of melt bottoms. The chemical and mechanical properties can only be guaranteed for new material. Melting with open flame (acetylene / oxygen) and inductive melting: Once the cylinders are melted and the cast shadow falls across the molten metal, before the oxide skin begins to split, start the casting. **Maximum temperature for casting: 1450 °C.**

**Devesting and cleaning**

Let the muffle cool down to room temperature (approx. 20 °C), do not quench with water. Put the cooled muffle into water to avoid dust generation during the devesting. Sandblast the surface with 110 µm of aluminium oxide with 3-4 bar, then clean with a steam cleaner.

**Soldering / Laser welding**

We recommend commercially available nickel based solders. Nicor parts should not be soldered with gold or palladium solder. Nicor is also ideally suitable for laser-welding (commercially nickel based laser wire).

**Preparation before ceramic veneering**

The frameworks can be elaborated with standard carbide cutters, look for smooth transitions and avoid overlapping material. Please use the same cutter for one alloy to avoid contamination. The minimum thickness of the prepared coping should not be less than 0.3 mm. It's recommended to sandblast the frames with minimum 110 µm of aluminium oxide with 3-4 bar and clean with steam cleaner. Oxide firing is not mandatory but can be done as an option for 5 minutes at 980 °C with vacuum (cleaning firing). The frame needs to be sandblasted with aluminium oxide 110 µm and 3-4 bar to remove the present oxide layer thoroughly. In the end the cleaning by steam cleaner is mandatory. If you use a ceramic bonder please consider the instruction for use of the manufacturer.

**Handling conditions / Safety**

**Metal dust is harmful to health. Use when grinding and sandblasting dust extraction and respirator with filter FFP3 – DIN EN 149.**

**This product contains Nickel! May produce an allergic reaction. Safety data sheet available on request. Only for professional user!**

**Contraindications and side effects**

If the instructions are observed during the production processes, incompatibilities with NiCr alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the contra-indications and side effects. Any serious incident that involve the product must be reported to the manufacturer and the competent authority in the accorded country.

**Single-use**

The chemical and mechanical properties can only be guaranteed for new material.

**Disposal Instructions**

Consult the material safety data sheets or national regulations for disposal. Dispose of Nicor residues and dust in an environmentally friendly manner. Grinding dust must not enter groundwater, water bodies or sewers. Address waste exchanges for recycling.

**Storage conditions**

Temperature, humidity or light has no effect on the product properties.

**Quantity**

Please consider the label on the package.

**Dimension**

Ø 8 mm x 15 mm

**Chemical composition (typical values)**

Ni %	Cr %	Mo %	Si %	Others %
61.4	25.9	11	1.5	< 0.1

Determination, extent and validity of the chemical composition according to DIN EN 10 204 – 3.1

**Typical technical data (After casting)**

Yield strength 0.2%	340 MPa
Elongation	26.5%
Tensile strength	550 MPa
E-module	178 GPa
Density	8,43 g/cm <sup>3</sup>
Corrosion resistance	< 200 µg/cm <sup>2</sup>
Tarnish resistance	yes
Hardness	185 HV 10/30
CTE (25-500 °C)	~ 14,1 x 10 <sup>-6</sup> K <sup>-1</sup>
Melting range (Solidus / Liquidus)	1325 °C / 1350 °C

**Applied Norm:**

DIN EN ISO 22674

*Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.*

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