


Tizian Blanks

- für Pantographen
- for pantographs
- pour pantographes
- per frese pantografe
- для пантографов
- für CAD/CAM
- for CAD/CAM
- pour CAD/CAM
- per CAD/CAM
- для CAD/CAM

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Processing instructions

Description

Tizian Blanks are composed of tetragonal zirconium dioxide stabilized with yttrium oxide. It can be obtained in different blank sizes and permits the production of full ceramic tooth replacements. The materials are partially sintered and optimised for use in copy milling or CAD/CAM procedures.

Applications

The blanks are suited to the production of

- Primary conus and telescope crowns
- Crown frameworks in front and side tooth areas
- Bridge frameworks in front and side tooth areas
(Anterior: max 3 pontics; Preoudor / molar: max 2 pontics)
- Abutments, inlays

Composition

Component
zirconium dioxide
yttrium oxide
hafnium oxide
aluminium oxide
silicon oxide

quantity [% weight]
< 96
> 4
> 1
< 1
< 0.02

Technical data

Colour	white
Density (pre-sintered)	approx. 3.1 g/cm ³
Impact green strength (pre-sintered)	approx. 55 MPa
Density (dense sintered)	approx. 6.09 g/cm ³
Breaking strength (dense sintered)	approx. 1200 MPa
E-modulus (dense sintered)	approx. 210 GPa
Vicker's hardness (dense sintered) HV10	approx. 1200
Open porosities	0%

Schematic representation of the processing of Tizian Blanks

1. Form grinding / Shape grinding

The process depends on the chosen to avoid procedure. Follow instructions on preparation!

2. Framework dimensions

Wall thickness must be at least 0.7 mm prior to sintering. Framework wall thickness should never fall below 0.4 mm after dense sintering. Following dense sintering, the connecting cross-sections should not fall below 9 mm² in the region of the anterior teeth and 12 mm² in the region of the side teeth.

It is important that when sintering several workpieces additional support is used to stabilise the pieces and avoid distortion.

3. Sintering in the high-temperature furnace

Prior to dense sintering, the frameworks must be freed from dust caused by grinding. The framework milled approx. 25% greater in size is shrunk to its original by dense sintering. A suitable furnace is required for this (to 1600 °C), e.g. Tizian Furnace.

Sintering program:

Start: room temperature
Heating rate: 8 °C/min.
Final temperature: 1480 °C
Dwell time: 120 min.
Cooling rate: 8 °C/min. (recommendation 6 °C/min.)
Opening the furnace at 300 °C

Work pieces spanning larger areas must be well supported. This is the only way of preventing distortion during sintering. The material obtains the desired physical properties through sinter firing, e.g. a maximum strength of approx. 1200 MPa.

4. Framework fitting

- Zirconium dioxide frameworks should only be processed under cooling with water!
- Use only low rpms, exert minimum pressure on the material and ensure cooling is sufficient (to avoid of thermal damage!).
- Use only clean, news diamond-based instruments.
- **We do not recommend the use of blasting media for surface conditioning.** In case you do not wish to dispense with the use of blasting media, then do not exceed a working pressure of max. 2 bar (aluminium oxide max. 110 µm).

5. Veneering

Veneering should only be carried out using a veneering ceramic suitable for use on zirconium dioxide, e. g. Tizian. Zirconium dioxide frameworks must be fully veneer. The only exceptions are scalloped line or protective back plates for stabilization. These must be prevented from coming into contact with the antagonist.

6. Fixing

These restorations can be fixed in place both conventionally (glass ionomer or zinc phosphate cement) as well as adhesively. Provisional fixing in place is possible but should only take place for a limited time.

Note

The open processing of partially sintered zirconium dioxide could lead to alveolar dust. A protective mask must therefore be worn during processing.

Storage

The material should be stored in a dry place at 10 - 25 °C.

Waste disposal

The material can be disposed of in household waste.

Guarantee

The technical instructions on processing, are based solely on our own experience and can therefore only be regarded as guidelines no matter whether given written, spoken or by practical demonstrations.

Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-indications / interactions

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

Note

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

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Schütz Dental GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany