# Tizian Blank Modell

Schütz Dental GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany  
Telefon +49 (0) 6003 814-0 · Telefax +49 (0) 6003 814-906  
info@schuetz-dental.de · www.schuetz-dental.de

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## Instructions

### Description / Indications

Tizian Blank Modell are indicated for the production of dental models. The blanks are made of a special, highly compressed ureol with exceptional edge strength and abrasion resistance. Tizian Blank Modell exhibit an excellent machinability. They are suited to the production of model stumps and model bases.

### Important information

Tizian Blank Modell are not permitted for patient cases!

### Physical properties

chemical composition	ureol (reaction product on poly-urethane basis)
color	beige
Shore D-hardness	75–85
density	1.0 g/cm <sup>3</sup>
coefficient of thermal expansion	50–55* 10 <sup>-6</sup> K <sup>-1</sup>
dimensional stability under heat (HDT)	70–80 °C
compressive strength	45–50 MPa
flexural strength	50–55 MPa

### Processing

The optimal ambient temperature for milling Tizian Blank Modell is 20–25 °C. Please always follow the instructions supplied by the manufacturer of your milling machine.

The milling parameters supplied below are to be regarded as guidelines only. The effective optimum values depend on the milling machine, the CAM software as well as on

the milling geometries used. If uncertain, please use the same parameters that you use for PMMA blanks.

### Please note!

When milling Tizian Blank Modell, respirable dust can be discharged which may damage the lungs or lead to irritations of the eyes and skin. Therefore, the material must only be processed while using a proper and efficient extractor system.

Please follow the below listed safety advice:

- Always process Tizian Blank Modell underneath a local extractor system.
- Ensure that the material does not get into your eyes or come into contact with your mucous membranes.
- Do not inhale milling dust.
- We recommend to wear protective masks and goggles.

### Storage

Store in a dry place at temperatures between 10 °C and 25 °C (50 °F and 77 °F). Protect from direct sunlight and from humidity.

### Disposal

The material can be disposed of with normal household trash.

### Guarantee

Our processing instructions, regardless of whether they are issued verbally, in writing or in the form of a practical demonstration, are based upon our own experience and can therefore be regarded solely as guidelines.

parameters	roughening	coarse finishing	fine finishing
recommended mill diameter (mm)	2.5–4.0	1.5–2.5	0.6–1.0
feed rate (mm/s)	60	40	80
feed depth (mm)	max. 2.0	–	–
path distance (mm)	–	0.3–0.5	0.1–0.2
milling allowance (mm)	0.1	0	0

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