# dialog Vario

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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## Instructions for use

dialog Vario is a light curing composite for the crown and bright technology.

### Range of applications

#### **Fixed restorations**

- Crowns, bridges, adhesive bridgework
- Inlays, onlays, implant supra-structures, laminate veneers, longterm temporary restorations

#### **Removable restorations**

- Conical and telescopic crowns
- Attachment matrices
- Customizing acrylic denture teeth

### Composition dialog Vario

**Anorganic filler** (74 % by weight and 66 % by volume)

Glass filler (average particle size: 0.7 µm)

Pyrogenic silica (average particle size: 0.04 µm)

**Monomers** (25 % by weight)

Urethane dimethacrylate

1,4-butanediol dimethacrylate

BisGMA

**Additives** (1 % by weight)

Initiators, stabilizers, pigments

### Composition dialog Vario Flow

**Anorganic filler** (57 % by weight, 51 % by volume)

Glass filler (average particle size: 0.7 µm)

Pyrogenic silica (average particle size: 0.04 µm)

**Monomers** (43 % by weight)

Urethane dimethacrylate

1,4-butanediol dimethacrylate

**Additives** (1 % by weight)

Initiators, stabilizers, pigments

### Directions for use

The crown or bridge framework is waxed up using standard techniques and provided with retentions.

## Sebond Smart

### **Application**

Sebond Smart is used as a bonder between acrylics / composites and all types of metal.

Following processing, the area being veneered is sand-blasted with aluminium oxide (50-110 µm, 2 bar). The area thus activated is then cleaned with oil-free compressed air. Sebond Smart is applied

generously with a brush and allowed to flash off for approx. 1 min. A second layer can be applied if necessary. Do not blow off Sebond Smart with compressed air!

### Storage

Sebond should be stored at 10-25 °C. Close the bottle properly immediately after use and protect against the effects of direct light.

### Shelf life

The maximum shelf life is printed on the label of each bottle. Do not use after the expiry date.  
Each bottle contains: 5 ml

## dialog Vario paste opaquer

The paste opaquer is initially applied like a wash using a stiff, short-haired brush. Apply the opaque in two to three thin coats to achieve full coverage. Polymerize each coat separately. An opaque coat without full coverage will have a negative influence on the color effect. The viscosity of the opaque can be optimized by stirring on a mixing block. Particular attention should be paid to the application of very thin coats in retention areas.

For polymerization times, see table 2.

### Sources of error / error prevention

The surfaces of the framework must be cleaned of any grease and polish residue. The opaquer should be applied only in very thin layers.

Opaque coats that are too thick prevent optimum polymerization and thus weaken the bond. Extension of the polymerization time does not result in a greater curing depth.

### Working time for the materials

Depending on the lighting conditions 1-3 min.

### Storage

The material should be stored at 10-25 °C. Close the syringe properly immediately after use and protect against the effects of direct light.

### Shelf life

The maximum shelf life is printed on the label of each syringe. Do not use after the expiry date.  
Each syringe contains: 3 g

## dialog Vario and dialog Vario Flow

### Dentin material

Extrude the dentin material from the syringe and apply it to the entire area to receive the facing, tapering a thin layer toward the incisal area.

Depending on the shade, use a spatula to contour a slightly corpulent dentin core.

Carry out an intermediate polymerization, see table 2.

If a dialog veneer is coated using thicknesses exceeding 2 mm, the intermediate polymerization time must be extended to 8 min.

### Modelling aid

To make the mass somewhat smoother and prevent it from sticking to the modelling tool, a drop of modelling aid can be added to the composite and mixed in with it.

### Incisal material / Incisal occlusal material

Extrude the desired incisal material from the syringe and apply it to the incisal half of the veneer, over the dentin layer. The incisal material is used to provide the veneer with its final contours. Once the veneer has been contoured as desired, polymerize it.

For polymerization times, see table 2.

The entire surface of the veneer should be coated with expers Gel prior to final polymerization. This prevents formation of a further smear layer and facilitates trimming and finishing.

### Fixation of surface properties

After polymerization as described previously, the veneer exhibits an increased yellow value due to the catalyst. In order to achieve the final shade and retain it, the restoration must be polymerized again.

### Finishing the faging

Silicone polishers (square-edge or knife-edge wheels or cylinders), tungsten carbide cutters and diamond instruments are suitable for finishing.

### Polishing

The dialog veneer is polished with a handpiece mounted goat hair brush, polishing paste (e.g. dialog Vario polish) and soft woolen buffs. Careful finishing and polishing of the surface is essential, if optimum results are to be achieved. What's more, it can prevent deposits to a great extent which could cause discoloration, such as nicotine, caffeine etc.

**Please note:** Like all acrylics, dialog produces fine dust while being trimmed. It is advisable to work over a dust extractor and to wear protective goggles.

### Correction and repair

Roughen the veneers up to 2 mm around the margin of the area to be corrected or repaired, brush Bonding Fluid onto it and light cure it to create a new smear layer.

For polymerization times, see table 2, „intermediate polymerization“.

Should the polymerized layer of Bonding Fluid show a whitish appearance, either the layer is too thick or the material has been polymerized excessively. In this case, it must be removed. Repeat the previous work-

step, but reduce the layer thickness and the polymerization time. Next the appropriate dialog material is applied and polymerized.

### Chipping of veneers

- Wax-up the framework properly to ensure that the dialog materials are always supported regardless of the prevailing occlusal conditions.
- Avoid precontacts.
- Apply the individual materials in the correct thickness and polymerize for the correct periods.
- Add retention beads.
- After letting Sebond Smart dry for 1 minute, immediately apply the opaquer.

### Entrapped air

- Rotate the plunger to extrude the paste out of the syringe and scrape it off – do not use an instrument to remove the paste from the syringe.
- Apply sufficient material for the veneer and spread it.
- Do not mix dialog materials, rather apply coats on top of each other.

### Discoloration and plaque deposits

- Check the luminous power of the light curing units.
- Polish the surface thoroughly to seal it.
- Ensure that the restorations are positioned correctly in the light curing units.
- Apply the dialog materials in the correct thickness and adhere to the polymerization times.

### Guarantee

Our technical instructions, regardless of whether they are provided verbally, in writing or during practical demonstrations, are based on our own experience and should only be considered guidelines.

As our products are subject to continued development, we reserve the right to modify them.

### Processing time of the materials

1-3 minutes, depending on the light conditions.

### Curing depth of polymerization

The curing depth of polymerization is 2.0 mm. To ensure optimum physical properties, no layer should be more than 2.0 mm thick. A light curing unit with an emission spectrum of 310-500 nm must be used. The requested physical characteristics can only be obtained with immaculate lamps. Therefore, a regular check of the light intensity according to the values quoted by the manufacturer is necessary.

### Storage

To ensure the product stability, store the material between 10 °C and 25 °C. If, during high temperature periods outside, dialog is kept in a refrigerator, it must be removed well before use to resume its original contouring properties. Close the syringe tightly after use and protect

it against direct light. Turn back spindle about one turn, in order to keep the material from leaking out.

### Shelf-life

The maximum shelf-life is printed on the label of each syringe. Do not use after the expiry date.

## dialog Vario Occlusal / Veneer materials

dialog Vario Occlusal / Veneer materials are light curing composites for occlusal veneers, palatal surfaces, for building up cusps and for full-coverage veneers.

Available are 8 dialog Occlusal incisal materials in the colors SR (incisal red), SL (incisal light), clear (CL) and blue (BI), S57, S58, S59 and S60.

### Composition

#### Monomer matrix

Diurethane dimethacrylate; Bis-GMA; 1,4-butanediol dimethacrylate

#### Total filler content

Inorganic filler: 75 % by weight (52 % by volume); glass filler (mean particle size: 0.7 µm); pyrogenic silicic acid (mean particle size: 0.04 µm)

### Indications

- Cosmetic adjustments
- Indirect laminate veneers
- Vestibular veneers on temporary full and partial dentures
- Posterior veneers
- Inlays
- Onlays
- Full-coverage veneers

### Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

### Contra-indications / interactions

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

### Note

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

## Interactions with other materials

As phenolic substances (such as eugenol) inhibit polymerization, do not use materials containing these substances.

### Dispensing and application

These occlusal materials must **not be applied directly** to metal or opaquer. Once the bonding agent and opaquer have been applied to the metal framework, it must be coated with a layer of dialog Vario as described in the dialog instructions. This flexible intermediate layer absorbs loads and stress. The shade is then selected from the enamels and the occlusal materials built up. Do not exceed a layer thickness of 2 mm. Once they have been built-up, dialog Vario materials only require intermediate polymerization before applying dialog Occlusal.

## Indirect inlays (indirect method)

### Fabricating inlays

Pour the impression with high strength stone in the laboratory. Once the model has set, remove the impression. Block out the undercuts and apply an oil-free separating agent (e.g. Futura Sep, Schütz Dental) to the model.

Build up the inlay layer-by-layer on the model, starting with the approx- imal and deep occlusal areas. Maximum thickness of each layer: 2 mm. Shade modifiers (e.g. dialog Vario Chroma Flow, Schütz Dental) may be applied between the dentin and enamel to create special effects. Each layer must be cured individually with a standard light curing unit (see table 2). The finished inlay is then released from the die and polymerized. The occlusal surface should be trimmed with fissure burs. The inlay is polished using goat hair brushes, polishing paste and soft woolen buffs in a handpiece. Careful trimming and polishing of the surface is essential for optimum results and virtually prevents deposits (nicotine, caffeine etc.) accumulating on the surface and affecting the shade.

**Please note:** As with all composites, dialog Occlusal produces fine dust during trimming. It is advisable to use a dust extractor.

### Placing inlays, onlays, laminate veneers

Remove the temporary restoration and clean the cavity. Place a rubber dam, clean and dry the prepared tooth surfaces. Apply gentle pressure on the restoration to ensure that it fits precisely. Once the tooth has been cleaned, it is etched (e.g. Capo Etch, Schütz Dental) and the bonding agent (e.g. Capo Bond, Schütz Dental) applied and cured with a commercially available curing unit (refer to the separate instructions). A luting composite (e.g. Alphaslink Cem, Schütz Dental) is then applied as described in its manufacturer's instructions. The inlay is placed in the cavity and gently pressed into position.

**Caution:** If the inlay is more than 2 mm thick, a dual-curing composite should be used. Excess composite is removed and the remaining composite cured with a commercially available light curing unit as described in its manufacturer's instructions. Remove excess composite with fine grit diamonds and diamond finishing strips. Check the occlusion and adjust where necessary.

## Polishing

Trimming and polishing are carried out with the finishing and polishing set. We recommend the use of our dialog Vario polishing paste which is matched to this system.

**Caution:** Avoid entrapping bubbles. Individual layers should not be thicker than 2 mm. Due to the influence of oxygen in the atmosphere, a thin non-cured film (smear layer) remains on the surface of every layer. This smear layer ensures that the layers bond together chemically and must not be touched or contaminated with moisture.

We recommend applying expers Gel (Schütz Dental) prior to final curing. This gel prevents the material from reacting with oxygen and allows the surface to cure fully. The working time beneath a dental light is approximately three minutes, depending on the lighting conditions. If extensive contouring is required, the lighting intensity should be reduced or the material covered with an opaque foil.

## Curing

The intermediate light curing time for all shades is approximately 60-90 seconds per layer when using a commercially available light curing unit. The final shade and hardness is not attained until the final layer has been applied and cured (refer to table 1).

### Curing details

A light curing unit with an emission spectrum of 310-500 nm must be used for curing. As the required physical properties are only attained if the light curing unit is functioning properly, its luminous intensity must be checked regularly as described by the manufacturer.

### Usage and storage

Store at 10°C-25°C. Do not use after the expiry date printed on the syringe label. Ensure that the material is at room temperature before use. For use by dental professionals only. After dispensing material, withdraw the plunger slightly to prevent the nozzle becoming blocked. After dispensing material, close the syringe immediately and always keep it closed. If the material is not cured fully, it may discolor, its full mechanical properties will not be achieved and irritation of the pulp cannot be ruled out. This product was developed specifically for the range of indications described. It must be used as described in the instructions. The manufacturer cannot be held responsible for damage due to incorrect handling or processing.

Table 1

Color according to Vita**	Paste opaquer	Dentin	Incisal
A1	light	A1/B1	S 58
A2	light	A2/B2	S 58
A3	light	A3	S 59
A3,5	medium	A3,5	S 59
A4	medium	A4	S 60
B1	light	A1/B1	S 57
B2	medium	A2/B2	S 59
B3	light	B3/B4	S 59
B4	medium	B3/B4	S 59
C1	medium	C1/D2/D3	S 59
C2	medium	C2/C3	S 59
C3	dark	C2/C3/C4	S 59

C4	dark	C2/C3/C4	S 60
D2	light	C1/D2/D3	S 60
D3	dark	D2/D3	S 59
D4	dark	D2/D3/D4	S 60

\* Vita is a registered trademark of Vita Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen/ Germany.

### Recommendation for layering

- B2: On the edges, draw incisal all the way down to the shoulder area, but cover the shoulder area only with a relatively thin layer of incisal material.  
 C1: Completely cover tooth with incisal material.  
 C4: Modify dentin C3 with khaki/blue mixture. On the edges, narrow out the incisal material thinly towards the shoulder, apply a very thin layer to the tooth body.  
 D4: Modify dentin D3 with khaki/blue mixture. Only on the edges, draw the incisal material all the way down to the shoulder area, thinly cover tooth body with incisal material.

**Table 2**

Light curing device	Paste opaquer	dialog Vario / dialog Vario Flow			dialog Vario Occlusal		
		Intermediate polymerization	Final polymerization	Surface treatment	Intermediate polymerization	Final polymerization	Surface treatment
Spektra LED	1 min.	30 sec.	3 min.	none	30 sec.	3 min.	none
Spektra 2000 fast	3 min.	90 sec.	6 min.	5 min.	90 sec.	6 min.	5 min.
Spektra 2000	3 min.	90 sec.	9 min.	7 min.	90 sec.	9 min.	7 min.
Uni-XS	90 sec.	90 sec.	180 sec.	180 sec.	90 sec.	180 sec.	180 sec.
HiLite	90 sec.	90 sec.	180 sec.	180 sec.	90 sec.	180 sec.	180 sec.
Spektramat	2 min.	1 min.	5 min.	5 min.	1 min.	10 min.	10 min.
Labolight LV-II/III	1 min.	0,5 min.	9 min.	9 min.	0,5 min.	9 min.	9 min.
Solidilite EX / V	1 min.	90 sec.	9 min.	7 min.	90 sec.	9 min.	7 min.

**Table 3 – Paste opaquer**

Trouble	Cause	Remedy
material does not cure	layer too thick insufficient polymerization	only apply wash-like respectively in very thin layers Pobserve polymerization times, incorrect curing times/check curing light /exchange curing light, if necessary / curing light defective / exchange light source
chipping on opaquer /veneer	no retention opaquer not polymerized sufficiently /opaquer layer too thick	use retentions only apply wash-like respectively in very thin layers /incorrect curing times / curing light defective / exchange light source

**Table 4 – dialog Vario /dialog Vario Flow / Occlusal**

Trouble	Cause	Remedy
material does not cure	layer too thick insufficient polymerization	maximum layer thickness 2 mm observe polymerization times /check curing light /exchange light source, if necessary
chipping on veneer	opaque not polymerized sufficiently / opaquer layer too thick	only apply wash-like respectively in very thin layers /incorrect curing times / check curing light /exchange light source, if necessary
material inside the syringe too hard	syringe was stored at temperatures below 10 °C	store material at room temperature or let the material reach room temperature in time for the treatment
smear surface	insufficient polymerization	observe polymerization times /check curing light /regular maintenance
shade does not match	incorrect layering	observe color table / too much or too little incisal material was used
material is impossible to polish	insufficient polymerization dull surface	observe polymerization times /check curing light /exchange light source, if necessary insufficient polishing / use appropriate polishing paste
shade of the veneer appears too yellow	insufficient polishing incorrect layering	observe polymerization times /check curing light /exchange light source, if necessary observe color table

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Schütz Dental GmbH