# **Tizian Blank NEM**

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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## Instruction for use

### **Intended use**

Dental alloys and metals are used for the fabrication of fixed, removable or combined dental restorations.

### **Indication**

Crowns and bridges up to 16 units (max. 4 pontics in the anterior and max. 3 pontics in the posterior region), implant-supported superstructures, bars.

### **Contraindication**

- All indications not listed under Indication.
- In case of known allergic reactions to any of the ingredients.

### **Milling**

Tizian Blank NEM is intended for CNC milling machines. Please follow the instructions and parameters of the respective manufacturer of CAM Software and the CNC milling machine.

### **Framework design**

Design the objects using suitable CAD software, taking into account the rules of dental technology. Please consider an anatomically reduced framework design for the veneering with ceramic. The wall thickness should not be less than 0.3 mm. Ensure that the connector cross-section (6-9 mm<sup>2</sup>) is adequate for bridge links. Sharp edges and undercuts should be avoided.

### **Cutting out the frameworks from the blank**

Remove the milled frameworks with suitable hard metal mills or cutting discs and grind down the supports.

## **Soldering / Laser welding**

Tizian Blank NEM can be used with all suitable solders for CoCr alloys. Tizian Blank NEM parts must not be soldered with gold or palladium solder. Tizian Blank NEM is also ideally suitable for laser-welding.

## **Preparation of the surface before ceramic veneering**

The frameworks can be elaborated with standard carbide cutters, look for smooth transitions and avoid overlapping material. Please use the same cutter for one alloy to avoid contamination. The minimum thickness of the prepared coping should not be less than 0.3 mm. It's recommended to sandblast the frames with minimum 110 µm of Aluminium oxide with 3 - 4 bar and clean with steam cleaner. Oxide firing is not mandatory but can be done as an option for 5 minutes at 980 °C (cleaning firing). The frame needs to be sandblasted again to remove the present oxide layer thoroughly. Subsequently, clean the framework with a steam jet. If you use a ceramic bonder please consider the instruction for use of the manufacturer.

## **Handling conditions / Safety instructions**

Metal dust is harmful to health. When grinding and sandblasting use dust extraction and respirator with filter FFP3 – DIN EN 149.

## **Contraindications and side effects**

If the instructions are observed during the production processes, incompatibilities with non-precious dental alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the contra-indications and side effects.

## **Disinfection of the dental prosthesis before insertion**

Workpieces from the dental laboratory must be subjected to immersion or spray disinfection and then rinsed under running water before insertion into the patient's oral cavity.

## **Single-use**

Once milled, blanks must not be processed for further use as a medical device.

## **Disposal instructions**

Please dispose of metal residues and dust in an environmentally friendly manner. Do not allow waste to enter groundwater, water or sewage systems. Contact waste exchanges for recycling. Outer packaging can be disposed of in paper waste.

## **Storage conditions**

Temperature, humidity or light have no effect on the product properties.

## **Quantity**

Please consider the label on the package.

## **Description**

Cobalt-based non-precious metal dental milling ceramic alloy, type 4.

## **Dimension**

Ø 98.3 mm: 8 mm, 10 mm, 12 mm, 13.5 mm, 15 mm, 16 mm, 18 mm, 20 mm

## **Chemical Composition** (typical values)

<b>Co %</b>	<b>Cr %</b>	<b>W %</b>	<b>Si %</b>	<b>Mn %</b>	<b>Fe %</b>
<b>61.65</b>	<b>27.75</b>	<b>8.45</b>	<b>1.61</b>	<b>0.25</b>	<b>0.2</b>

## Typical technical Data

Yield strength 0.2%	440 MPa
Elongation at break	14%
Tensile strength	819 MPa
E-module	210 GPa
Density	8.55 g/m <sup>3</sup>
Corrosion resistance	< 200 µg/cm <sup>2</sup>
Hardness*	275 HV 10/30
CTE (25 - 500 °C)	14.1 x 10 <sup>-6</sup> K <sup>-1</sup>
Max. firing temp.	~ 980 °C

\*with appropriate sample preparation

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Schütz Dental GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany