# FuturaBasic Cold

---

CE 0297

Schütz Dental GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany  
Tel.: +49 (0) 6003 814-0 · Fax: +49 (0) 6003 814-906  
info@schuetz-dental.de · www.schuetz-dental.de

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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## Instructions for use

### Intended use:

Materials for the fabrication of partial dentures, full dentures, framework-supported dentures using different techniques such as the injection technique, the flow technique and the flask technique.

### General product description:

Cold-curing acrylic on methyl methacrylate basis.

### Composition:

- Powder: pearl polymer made of polymethyl methacrylate (PMMA), pigments, initiators
- Liquid: methyl methacrylate, tetramethylene dimethacrylate, initiators, stabilizers
- Futura Basic Cold is free of cadmium.

### Indications:

- For the production of partial and full dentures in the flow technique.
- For repairs, extensions and relinings of existing restorations in the flow technique.
- For the production of splints
  - occlusal splints
  - reflex/anti-grinding splints

### Contraindications / interactions with other agents:

Do not use in case of a known allergy against one of the components.

### Application / processing:

#### Preparatory steps:

Produce plaster model and wax-up as usual.

The prefabricated teeth set up in the wax are fixed with plaster or silicone (stent/flask). Subsequently, the wax is eliminated as usual.

Roughen the acrylic teeth, provide them with undercuts and coat them with monomer or HS Cross-Liquid (Schütz Dental) to ensure that they are securely retained. The surfaces of the plaster models are insulated with the alginate separating agent FuturaSep Plus (Schütz Dental) according to the manufacturer's instructions.

#### Mixing ratio:

10 weight parts powder, 7 weight parts liquid. Veins can be added, if desired. Please note the separate manufacturer's information and instructions.

Approx. **20 g powder** and **14 g liquid** are required for a full denture.

### Production of partial and full dentures as well as of implant-supported restorations and framework completions in the flow technique:

Mix powder and liquid (according to the mixing ratio). After allowing a doughing time of 15 sek., stir Futura Basic Cold with a spatula to produce a homogeneous mixture. Fill the die as usual.

After mixing at room temperature (21 °C / 69 °F), Futura Basic Cold is pourable for ca. 2 min. After 3 min., the material enters the plastic phase.

After the plastic phase has begun, the resin no longer flows out of the stent and can be modeled.

After 5 min. (at a room temperature of 21 °C / 69 °F), at the very latest, the restoration must be placed in the pressure pot.

Polymerize Futura Basic Cold in a pressure polymerization system for 15 min. at a pressure of 2–5 bar and a water temperature of 45 °C / 113 °F.

### Repairs and relining of existing acrylic appliances in the flow technique:

If necessary, produce a mold or a stent.

Clean the acrylic surfaces to be repaired, roughen them and smooth out the edges. Wet with monomer or with HS Cross-Liquid (Schütz Dental).

Mix powder and liquid (according to the mixing ratio). After allowing a doughing time of 15 sek., stir Futura Basic Cold with a spatula to produce a homogeneous mixture. Let the acrylic flow in.

Polymerize Futura Basic Cold in a pressure polymerization system for 15 min. at a pressure of 2–5 bar and a water temperature of 45 °C / 113 °F.

### Finishing the material:

After polymerization, the restoration is lifted from the model and it is finished with suitable tools.

Silicone polishers (wheel, lentil, roll) as well as hard metal cutters or diamond grinding tools are suited for finishing the material.

Prepolish with sandpaper or pumice.

**Polishing:**

Polish the material with goat hair brushes, polishing paste and soft woolen wheels. A careful surface treatment and polish are a precondition for an optimal result.

**Production of splints:****Preparatory steps:**

Production of a plaster model.

Blocking undercuts with suitable materials.

Carve splint in wax as desired.

Produce impression, fit it with sprue canals and remove wax.

Insulate plaster model with suitable materials (e. g. FuturaSep, Schütz Dental).

**Processing:**

Mix powder and liquid (according to the mixing ratio). After allowing a doughing time of 15 sek., stir Futura Basic Cold with a spatula to produce a homogeneous mixture.

Pour acrylic dough into mold as usual.

These work steps must be completed after 6 min. at the very latest.

Polymerize Futura Basic Cold in a pressure polymerization system for 15 min. at a pressure of 2–5 bar and a water temperature of 45 °C / 113 °F.

**Finishing the material:**

After polymerization, the splint is lifted from the model and it is finished with suitable tools.

Silicone polishers (wheel, lenticle, roll) as well as hard metal cutters or diamond grinding tools are suited for finishing the material.

Prepolish with sandpaper or pumice.

**Polishing:**

Polish the material with goat hair brushes, polishing paste and soft woolen wheels. A careful surface treatment and polish are a precondition for an optimal result.

**Repairs:**

Produce a die or a stent, if necessary.

Clean the acrylic surfaces to be repaired, roughen them and smooth out the edges. Wet with monomer or with HS Cross-Liquid (Schütz Dental). Mix acrylic according to the mixing ratio and let the acrylic flow in.

Polymerize Futura Basic Cold in a pressure polymerization system for 15 min. at a pressure of 2–5 bar and a water temperature of 45 °C / 113 °F.

**Troubleshooting / FAQ list:**

| <b>Problem</b>                                | <b>Cause</b>                                                                                                               | <b>Remedy</b>                                                                                                                                                             |
|-----------------------------------------------|----------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| acrylic does not cure properly                | incorrect mixing ratio                                                                                                     | <ul style="list-style-type: none"><li>– observe instructions for use</li><li>– observe mixing ratio</li></ul>                                                             |
| acrylic shows a whitish tint along the bottom | <ul style="list-style-type: none"><li>– incorrect or improper insulation</li><li>– model was not soaked in water</li></ul> | <ul style="list-style-type: none"><li>– insulate properly</li><li>– check insulation agent, if necessary</li><li>– soak model in lukewarm water for ca. 10 min.</li></ul> |
| acrylic forms bubbles along the bottom        | <ul style="list-style-type: none"><li>– model was not soaked in water</li></ul>                                            | soak model in lukewarm water for ca. 10 min.                                                                                                                              |
| acrylic surface white /spongy                 | too much time passed until the acrylic was injected or until the acrylic was put in the pressure pot                       | observe processing times                                                                                                                                                  |
| color deviations of repairs                   | incorrect mixing ratio                                                                                                     | observe mixing ratio                                                                                                                                                      |
| acrylic fractures when lifted                 | improperly insulated                                                                                                       | check insulation                                                                                                                                                          |
| acrylic/prosthesis shows whitish streaks      | <ul style="list-style-type: none"><li>– improperly mixed</li><li>– powder content too high</li></ul>                       | <ul style="list-style-type: none"><li>– mix acrylic dough thoroughly with a spatula</li><li>– observe mixing ratio</li></ul>                                              |
| acrylic polymerizes too quickly               | <ul style="list-style-type: none"><li>– temperatures too high</li><li>– work steps too long</li></ul>                      | <ul style="list-style-type: none"><li>– observe temperatures</li><li>– observe working range</li></ul>                                                                    |
| splint fractures                              | design too thin /unsuitable                                                                                                | ensure a sufficient layer thickness (min. 2–3 mm)                                                                                                                         |

**General notes:**

For deviations from 21 °C room temperature: Higher temperatures reduce the above times, lower temperatures increase them.

To further reduce the risk of intolerance reactions, we recommend storing the orthodontic appliance in lukewarm water for at least 12 hours before insertion.

**Storage:**

Storage temperature at 10–25 °C / 50–77 °F.  
Close container carefully.

**Shelf-life:**

The maximum shelf-life is printed on the label of each packaging. Do not use after the expiration date.

**Side-effects:**

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (e. g. allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

**Contra-indications / interactions:**

If a patient has known hyper-sensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, the composition of the medical device we supply is available upon request. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth before using the product.

**Note:**

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

**Disposal:**

Powder and liquid must be disposed of as hazardous waste in accordance with official regulations.

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