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
If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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List of contents

Page

IMPLA – Implant System Micro Retention	3–6
IMPLA – Implant System Cylindrical	7–10
IMPLA – Mini-Implants conetop/balltop	11–12
General notes on implantation	13–14
IMPLA locking screws (covering and healing caps) and gingiva formers	15–16
IMPLA Abutments	17–18
Locator® Implant Abutment System	19–22
HSL Abutment Connectors	23–24
Multi-Unit Abutments	25–26
IMPLA PreFace Abutments	27–28
IMPLA Aesthura Abutments	29
Disinfection, cleaning and sterilization of surgical and prosthetic tools and instruments of the IMPLA system	30–33
Application and safety instructions	34
Manufacturer’s information on rotating surgical instruments Standardized manual reprocessing of rotating instruments	35–36
Information on storage, disinfection, cleaning procedures and sterilization of rotating instruments	37

IMPLA – Implant System Micro Retention

IMPLA implants are made of pure titanium (ISO 5832-2). The IMPLA system comprises surgical, prosthetic and laboratory components and instruments.

Important notes, warning

IMPLA implants are part of a complete system and may only be used with the original components and instruments of the system. The use of components which are not part of this system voids the guarantee and rules out parts being replaced. The IMPLA system may only be used by dentists, doctors and dental technicians familiar with dental implantology including diagnosis, preoperative planning, surgical techniques and prosthodontics.

Consultation is done in written form or by electronic media as well as with demonstrations during special workshops. All information is based on the state of the art and latest scientific findings known at the time of launching the product. However, we urge the users of our products to check the suitability of the product in each case. This is especially important for techniques not expressly recommended by the manufacturer. In case of doubt, we encourage our customers to confer with Schütz Dental Group/Section IMPLA.

Reuse and sterilization

All abutments and components are for single use only and must not be reprocessed. Any reuse involves the risk of cross-contamination! All instruments, abutments and components are packaged non-sterile and must be sterilized (EN 13060) prior to placement (see page 30–37). Using unsterile components increases the infection risk for the patient.

Indications

IMPLA implants Micro Retention are indicated for:

- restoring individual tooth gaps
- bridgework spanning interdental and free-end gaps according to general and special principles of prosthetic treatment guidelines.
- restoration of edentulous mandibles and maxillae

Further indications:

Implants with a diameter of 3.3 mm

These implants are an alternative in case of limited jaw ridge width. They are, because of their lower mechanical strength when compared with the standard diameter (4.2 mm), intended for the following indications:

- a) restoring individual anterior teeth of the mandible
- b) restoring lateral upper incisors
- c) four interforaminal implants to restore the edentulous mandible with a bar construction
- d) **Do not** use with ball anchors, magnets or telescopes

Implants with a diameter of 4.2 mm

The IMPLA standard implants with a diameter of 4.2 mm are used for all indications of oral enossal implants for mandibles and maxillae. A prerequisite is an adequate minimal amount of bone base with a width of about 5.5 mm.

Implants with a diameter of 5.3 mm

These implants are indicated for the use in the posterior tooth areas of mandible and maxilla (especially for the restoration of molars). The jaw ridge width should not be less than 6.5 mm.

The implants are primarily indicated for delayed implantations.

Contraindications

When selecting patients, the general contraindications for dental surgery must be taken into account. These include: Vascular diseases, uncontrolled diabetes mellitus and other metabolic or systemic diseases affecting bone regeneration, blood coagulation disorder, anticoagulative therapy, bone metabolism diseases.

Temporary contraindications: Acute and chronic inflammations, inadequate soft tissue present, chemotherapy or radiotherapy.

Oral contraindications:

- patient not really prepared to undergo total oral rehabilitation
- inadequate oral hygiene
- untreated parafunctions, e. g. bruxism
- inadequate amount of bone available
- inadequate interocclusal distance
- poor occlusion and/or articulation

Diagnostics / Patient information:

Diagnostics consist of a medical history, clinical examination, X-ray examination using OPG and, if necessary, CT. A medical check up is recommended. It is also advisable to perform biomechanical treatment planning to avoid overloading the bone, the implant or the prosthetic construction, taking into account the individual masticatory forces and occlusion. Implantations are elective procedures that require extensive patient education. Treatment alternatives and consequences of the implantations must be pointed out. In particular, the importance of oral hygiene and the costs of follow-up treatment must be pointed out. With regard to the type and scope of the information, reference is made to the relevant case law.

Side effects:

- Horizontal and vertical bone loss may come about even with a medically appropriate surgical procedure. The extent and the nature of bone loss cannot be estimated in advance.
 - Reversible or irreversible impairments or interactions are possible if the implants are placed near special structures (e.g. adjacent teeth, maxillary sinus or nerves).
 - In connection with dental implants, chronic pain may also occur in rare cases.
- To minimize any risk of infection, antibiotic prophylaxis may be recommended. In rare cases, allergies to the prosthetic components (immediate implantation) can also prevent successful implantation.

Important notes

When working with our implants and instruments intraorally, always use a protection against aspiration! Primary stability directly after the insertion of the implant is imperative for a successful healing. During the healing phase, IMPLA implants are generally not functionally loaded.

Surgical procedure

To create optimal conditions for the implant to heal, it is essential to treat the hard and soft tissues gently. It is advisable to use the countersink and bone thread cutter in order to ensure that the implant can be placed atraumatically applying as little torque as possible (maximum 30 Ncm). As thermal trauma may impede or prevent the implant healing, heat generation should be kept to a minimum as follows:

- recommended drill speeds:
2.0 mm Ø: approx. 800 rpm, 2.5 mm Ø: approx. 700 rpm, 3.15 mm Ø: approx. 500 rpm, Ø 4.2 mm, Ø 4.8 mm, Ø 5.3 mm: approx. 400 rpm
- recommended speeds: pilot drill Ø 1.2 mm and Ø 2.0 mm - approx. 800 rpm, enlargement drill Ø 3.15 mm - approx. 500 rpm, cavity drills of all Ø - approx. 400 rpm

- intermittent drilling and cutting using sharp drills and cutters (use 10 times only)
- cool the instrument externally with sterile physiological saline (NaCl) or Ringer's solution.

Important

To determine the optimum implant length and prevent nerve injuries, check in advance that adequate space is available (X-ray).

Placing the implant after preparing the implant bone supporting area

Removing the insertion aid

In order to remove the insertion aid (for implant 3.3 mm brown; 4.2 mm red; 5.3 mm blue), pull off the plastic implant carrier and loosen the vertical screw (blue) with the hexagonal screwdriver SW 1.2 mm and dispose of the vertical screw. Secure the implant against rotation with help of a counter wrench in cases of low primary stability and concomitant augmentation in the according alveole. Afterwards, the insertion aid can be lifted from the implant easily.

Closing the soft tissue

The cover / healing screw must be screwed onto the implant before the wound is closed with atraumatic suturing material.

Recovering the implant

Once the individual healing period (minimum of three months for the mandible and six months for the maxilla) has elapsed, the implant should be recovered and the cover / healing screw removed. In case of augmentation techniques, the healing period must be prolonged accordingly. The interior of the implant is then cleaned and a gingiva former suitable to the height of the mucous membrane is screwed in. Make sure that the gingiva former fits properly! The restoration should not be fitted to the implant until the soft tissue has healed fully and shows no irritations.

Impression taking

Open impression technique:

After the gingiva former has been removed, a forming aid with the appropriate diameter is fixed to the implant with help of a fixation screw. When taking the impression, the impression material is first syringed around the impression post. Once the impression material has cured, loosen the screw, remove the impression from the mouth together with the impression post and close the implant again with the gingiva former. Carefully screw an implant analog into the impression post and cast the model – preferably including a gingival mask. Once the cast has set, remove the impression post, respectively transfer coping, and replace it with the required abutment.

Closed technique:

Impression posts of the appropriate diameter are screwed onto the implant using a laboratory screw. Impressions are taken using a standard impression tray with the aid of transfer caps which are clicked onto the impression posts. The transfer caps remain in the impression. To produce a plaster model, they are fitted with a lab analog after the impression posts have been reinserted.

Laboratory components

When selecting the abutments, the angulation of the implant and the thickness of the gingiva must be taken into account. Those areas of the abutment which are in contact with the implant may not be blasted with abrasives or treated in any other manner.

Documentation

Each IMPLA implant can be traced via its lot number which is noted on a self-adhesive label for internal user documentation. Possible claims can only be asserted by providing this data.

Packaging and sterility

IMPLA implants are packed sterilely in a special clear acrylic container and are color-coded.

The chosen implant is removed from the sterile container with the insertion aid only directly before insertion and is inserted into the prepared bone cavity without touching the implant.

The covering and healing screw is attached to the back of the implant carrier. It can easily be removed from there with the hex screwdriver SW 1.2 mm.

Important reference data on implant length, implant diameter, article number, batch, CE certification, etc. are printed on the outer packaging made of cardboard.

General information

- All implants are for single use only.
- Each is packed sterilely in a closed container, please note the expiration date of sterility.
- The implant is sterile (radio sterilisation) as long as the packing remains unopened and undamaged. A subsequent sterilization by the user is not permitted by the manufacturer. A chemical sterilization of the implant (chemiclave) may not be carried out under any circumstances, as this could damage the implant surface. In case of damaged packaging, the implant is no longer sterile and must not be used anymore.
- Do not touch the implant.
- Store the implant in a dry place at room temperature. Avoid direct light incidence.
- System errors due to mixing up instruments and implants should be avoided. Please always check the color codings and labelings.
- Thermal damages of the jaw bone by rotating instruments (such as pilot drills, thread cutters or enlargement drills) must be avoided.
- It is advisable that the implantation is performed only after the patient information consultation has been completed, after the patient has signed the consent form and after the patient has undergone a preliminary medical examination for implantation.

Torque of screw connections

15 Ncm	horizontal screw bar screw cone adapter Multi Unit secondary prosthetic screw
20 Ncm	vertical screw cone adapter Multi Unit abutment 0° Multi Unit vertical screw Aesthura fixing screw
30 Ncm	Locator abutment ball connector

IMPLA – Implant System Cylindrical

IMPLA implants are made of medical grade 4 pure titanium (ISO 5832-2). The IMPLA system comprises surgical, prosthetic and laboratory components and instruments.

Important notes, warning

IMPLA implants are part of a complete system and may only be used with the original components and instruments of the system. The use of components which are not part of this system voids the guarantee and rules out parts being replaced. The IMPLA system may only be used by dentists, doctors and dental technicians familiar with dental implantology including diagnosis, preoperative planning, surgical techniques and prosthodontics.

Consultation is done in written form or by electronic media, with demonstrations during special workshops and by the implantology field staff. All information is based on the state of the art and latest scientific findings known at the time of launching the product. However, regardless of this, the user of our products is required to carry out a personal suitability test in the specific case and based on the indication. This is especially important for techniques not expressly recommended by the manufacturer. In case of doubt, we encourage our customers to confer with Schütz Dental GmbH/Section IMPLA.

Reuse and sterilization

All abutments and components are for single use only and must not be reprocessed. Any reuse involves the risk of cross-contamination! All instruments, abutments and components are packaged non-sterile and must be sterilized (autoclaved EN 13060) prior to placement (see page 30–37). Using unsterile components increases the infection risk for the patient.

Indications

IMPLA implants Cylindrical Implants are indicated for:

- restoring individual tooth gaps
- bridgework spanning interdental gaps and free-end gaps, in accordance with conventional and special principles of prosthetic treatment
- restoration of edentulous mandibles and maxillae

Further indications:

Implants with a diameter of 3.6 mm

These implants are an alternative in case of limited jaw ridge width. They are, because of their lower mechanical strength when compared with the standard diameter (4.5 mm), intended for the following indications:

- a) restoring individual anterior teeth of the mandible
- b) restoring individual lateral anterior teeth of the maxilla
- c) four interforaminal implants to restore the edentulous mandible with a bar construction

Implants with a diameter of 4.0 mm and 4.5 mm

The IMPLA standard implants with a diameter of 4.0 mm and 4.5 mm are used for all indications of oral enossal implants for mandibles and maxillae. A prerequisite is an adequate minimal amount of bone base with a width of about 5.3 mm and 5.8 mm.

Implants with a diameter of 5.5 mm

These implants are indicated for the use in the posterior tooth areas of mandible and maxilla (especially for the restoration of molars). The jaw ridge width should not be less than 6.8 mm.

The implants are primarily indicated for delayed implantations.

The implants with a length of 6.5 mm are generally not intended for single-tooth replacement or for accommodating a stand-alone prosthetic restoration.

Contraindications

When selecting patients, the general contraindications for dental surgery must be taken into account. These include: Vascular diseases, uncontrolled diabetes mellitus and other metabolic or systemic diseases affecting bone regeneration, blood coagulation disorder, anticoagulative therapy, bone metabolism diseases.

Temporary contraindications: Acute and chronic inflammations, inadequate soft tissue present, chemotherapy or radiotherapy.

Oral contraindications:

- patient not really prepared to undergo total oral rehabilitation
- inadequate oral hygiene
- untreated parafunctions, e. g. bruxism
- inadequate amount of bone available
- inadequate interocclusal distance
- inadequate occlusion and/or articulation

Diagnostics / Education:

Diagnostics consist of medical history, clinical examination, X-ray examination using OPG and CT, if necessary. It is recommended to have a medical check-up. It is also advisable to perform biomechanical treatment planning to avoid overloading the bone, implant or prosthetic construction, taking into account the individual occlusal forces and occlusion. Implants are elective procedures, about which extensive information must be provided to the patient. Treatment alternatives and consequences of the implantation must be pointed out. In particular, the importance of oral hygiene and the costs of follow-up treatment must be pointed out. With regard to the type and scope of the information, reference is made to the relevant case law.

Side effects:

- Horizontal and vertical bone resorption may occur even with a medically professional surgical procedure, and the extent and nature of bone resorption cannot be estimated in advance.
- Reversible or irreversible impairments or interactions are possible when implants are placed near special structures (e. g. adjacent teeth, maxillary sinus or nerves).
- Chronic pain may also occur in rare cases in connection with dental implants.

To minimize any risk of infection, antibiotic prophylaxis may be recommended. In rare cases, allergies to the prosthetic components (immediate implantation) can also prevent successful implantation.

Important notes

When using implants and instruments intraorally, it is generally important to ensure that they are secured against aspiration!

Primary stability directly after the insertion of the implant is absolutely necessary for a successful healing. During the healing phase, IMPLA implants are generally never exposed to a masticatory load.

Surgical procedure

To create optimum conditions for the implant to heal, it is essential to treat the hard and soft tissues gently. It is advisable to use the countersink in order to ensure that the implant can be placed atraumatically applying as little torque as possible (maximum 30 Ncm). The so-called „Dense Bone Drills“ have been designed especially for the cortical bone preparation. They are meant to reduce the insertion torque for longer implants. As thermal trauma may impede or prevent the implant healing, heat generation should be kept to a minimum as follows:

- recommended drill speeds:
2.0 mm Ø: approx. 800 rpm, 2.75 mm Ø: approx. 700 rpm, 3.0 mm Ø, 3.25 mm Ø and 3.4 mm Ø: approx. 500 rpm, 3.7 mm Ø, 3.9 mm Ø, 4.25 mm Ø, 4.9 mm Ø and 5.25 mm Ø: approx. 400 rpm
- intermittent drilling and cutting using sharp drills and cutters (use 10 times only)
- cool the instrument externally with sterile physiological saline (NaCl) or Ringer's solution.

Important

To determine the optimum implant length and prevent nerve injuries, check in advance that adequate bone is available (X-ray).

Placing the implant after preparing the implant bone area**Removing the insertion post**

In order to remove the insertion post (for implant 3.6 mm brown; 4.0 mm green; 4.5 mm red; 5.5 mm blue), pull off the plastic implant carrier and loosen the vertical screw (blue) with the hexagonal screwdriver 1.2 mm. Dispose of the vertical screw. In cases with concomitant augmentation or where the primary stability is low, secure the implant against rotation with help of a counter wrench. Afterwards, the insertion aid can be lifted from the implant easily.

Closing the soft tissue

The cover / healing screw must be screwed onto the implant before the wound is closed with atraumatic suturing material.

Exposure of the implant

Once the individual healing period (minimum of three months for the mandible and six months for the maxilla) has elapsed, the implant should be exposed and the cover / healing screw removed. In case of augmentation techniques, the healing period must be prolonged accordingly. The interior of the implant is then cleaned and the correct height of gingiva former for the oral mucosa is screwed in. Make sure that the gingiva former fits properly! The restoration should not be fitted to the implant until the soft tissue has healed fully and shows no irritations.

Impression taking**Open technique**

After removing the gingiva former, the impression aid of the appropriate diameter is screwed into the implant using a fixing screw. If an impression is subsequently taken with a customized impression tray, the impression aid is first coated with impression material. After the impression material has hardened and the fixing screw has been loosened, the impression together with the impression aid is removed and the implant is closed again with the gingiva former. A lab analog is carefully screwed to the impression post and then a model - preferably with a mucosa mask - is created. After the plaster has hardened, the impression aid or transfer abutment is removed and replaced by an appropriate prosthetic abutment.

Closed technique

Impression posts of the appropriate diameter are screwed onto the implant using a laboratory screw. Impressions are taken using a standard impression tray with the aid of transfer caps which are clicked onto the impression posts. The transfer caps remain in the impression and, after reinsertion of the impression posts, are fitted with lab analogs to produce a plaster model.

Laboratory components

When selecting the abutments, the angulation of the implant and the thickness of the gingiva must be taken into account. Those areas of the abutment which are in contact with the implant may not be blasted with abrasives or treated in any other manner.

Documentation

Each IMPLA implant is traceable by its lot number, which is noted on a self-adhesive label for internal user documentation. Any claims can only be asserted if this data is provided.

Packaging and sterility

IMPLA implants are packed sterile in a special clear container.

The chosen implant is removed from the sterile container with the insertion aid only directly before insertion and is inserted into the prepared bone cavity without touching the implant.

A special attachment on the bottom of the implant container holds the covering and healing screw which can be easily removed with the hex screwdriver SW 1.2 mm.

The outer cardboard box shows the user important data such as implant length, implant diameter, article no., lot and CE certification.

General information

- All implants are for single use only.
- Each implant is packed sterilely in a closed container, please note the expiration date of sterility.
- The implant is sterile (radiosterilization) as long as the packing remains unopened and undamaged. A subsequent sterilization by the user is not permitted by the manufacturer. A chemical sterilization of the implant (chemiclave) must not be carried out under any circumstances, as this could damage the implant surface. In case of damaged packaging, the implant is no longer sterile and may not be used anymore.
- Do not touch the implant.
- Store the implant in a dry place at room temperature. Avoid direct light incidence.
- System errors due to mixing up instruments and implants are to be avoided. Pay attention to the color codings and labelings.
- Thermal damages of the jaw bone by rotating instruments (such as pilot drills and extension drills) must be avoided.
- It is strongly recommended to perform implant surgery only after the patient consultation, after having the patient sign a consent form and after a medical check up. A surgery protocol is also recommended. Follow-up care with patient recalls to monitor the healing process and to support the patient's behavior after implantation is self-evident.

Tightening torque (torque) of screw connections!

10 Ncm horizontal screw connections

20 Ncm vertical screw connections

30 Ncm insertion of the implant, Locator Abutment screw connections

IMPLA – Mini-Implants conetop / balltop

Note:

Before performing implantology treatment, it is essential to be familiar with pertinent professional literature. It is absolutely necessary to have experience in surgery and in the assessment of the normal and pathological course of treatment.

The „General notes on implantation“ (see page 13–14) must be considered before treatment.

Description:

The one-piece IMPLA Mini-Implant is made of pure grade 4 titanium in accordance with ISO 5832-2.

Warning: If the implant breaks, it is possible that fragments remaining cannot be removed without damage to the surrounding bone tissue!

Implant threads: The IMPLA Mini-Implant has primary self-cutting threads, allowing primary stable insertion into the jaw bone while preserving bone material. The screw design and resulting compression stability together with the pilot drill and the extension drill ensure sufficient primary stability even in spongy bone tissue.

Implant surface: The implant was designed with a structured surface in accordance with the implant's application as a permanent implant (see indications).

Indications:

The IMPLA Mini Implant was developed for the following indications:

Mini implant „conetop“:

- Fixing full dentures
- Bar-supported dentures (on a minimum of 4 implants)
- Support implant for fixing bridges

Mini implant „balltop“:

- Fixing full dentures with help of an O-ring attachment

Contraindications:

Absolute:

- Patients whose systemic condition does not allow implantation.

Relative:

- Patients who are not capable of follow-up care for psychological or physical reasons and whose compliance is not anticipated
- Patients undergoing treatment for a malignant tumor
- Patients with uncontrolled diabetes
- Bruxism
- Heavy smokers
- Insufficiently available bone tissue and/or immediate vicinity to endangered structures (nerve, sinuses)

Diagnostics / Patient information:

Diagnostics consist of a medical history, clinical examination, X-ray examination using OPG and, if necessary, CT. A medical check up is recommended. It is also advisable to perform biomechanical treatment planning to avoid overloading the bone, the implant or the prosthetic construction, taking into account the individual masticatory forces and occlusion. Implantations are elective procedures that require extensive patient education. Treatment alternatives and consequences of the implantations must be pointed out. In particular, the importance of oral hygiene and the costs of follow-up treatment must be pointed out. With regard to the type and scope of the information, reference is made to the relevant case law.

Side effects:

- Even in a professional surgical procedure, potential horizontal and vertical bone loss must be taken into account as with all dental implants; the extent and kind of bone loss cannot be estimated in advance.
- Reversible or irreversible impairment or interaction is possible when implants are inserted in the vicinity of certain structures (e. g. adjacent teeth, maxillary sinus, or nerves).
- In rare cases, chronic pain can occur in association with dental implants.

To minimize the potential risk of infection, antibiotic prophylaxis can be recommended. In rare cases, allergies to the prosthetic components (immediate implantation) may prevent successful implantation.

Procedure:

The entire course of treatment has five stages:

1. Pilot drilling, 1.8 mm
2. Pilot drilling, 2.0 mm
3. Extension drilling, 2.5 mm
4. Initial implant insertion
5. Further insertion all the way to the implant's final position with help of the IMPLA screwdriver.

Point 1) Pilot drilling is performed at the optimal site (optionally also through the mucosa) with specification and control of the implant axis and length.

Points 2) and 3) The cavity is then finally prepared with the extension drill (diameter 2.5 mm, only for Mini conetop) according to the implant diameter 2.1 mm, 2.5 mm or 3 mm to be inserted. The distance to the definitive implants should not fall short of approx. 2 mm.

The treatment must be adapted to the bone situation:

D1–D2 => 1.8 mm => 2.0 mm => 2.5 mm drill in the lower jaw

D3–D4 and upper jaw => 1.8 mm => 2.0 mm => 2.5 mm drill with decreased drilling depth

Cooling with sterile saline solution is advisable.

Important: Drill alternately at a speed of 500–600 rpm.

The entire implant length should be predrilled with the 1.8 mm respectively 2.0 mm pilot drill.

Point 3) The preparation depth with the extension drill 2.5 mm should be adapted to the bone quality (see point 1).

Point 4) The packaging must not be opened until immediately before insertion of the implant. Remove the glass tube with the implant inside and remove the implant with the insertion aid. With help of the insertion aid, insert the implant into the cavity. Once an increasing torque is reached, the plastic insertion aid is removed.

Point 5) Use the IMPLA screwdriver to screw in the implant until the top edge of the structured implant area is level with the bone. If necessary, use the torque wrench to reach the final position.

Important: Because of the thread configuration of the implant, a torque of min. 25 to max. 50 Ncm should be reached, depending on the indication.

Notes:

- After implantation, the type of implant used and the LOT number must be recorded in the documentation (sticker).
- Implants must be stored in a closed container in a dry place. Before insertion, any contact with foreign substances must be avoided. The endosseous part of the implant must not be touched. In the event the implant, superstructure, or an accessory is swallowed, the location of the object must be determined and the appropriate medical steps are to be taken.
- Healing at the implant site must be monitored regularly, e. g. by x-ray. Implants that are infected or not osseointegrated must be removed promptly. The dentist must determine the time for removal.

Scope of delivery:

The IMPLA Mini-Implant is delivered sterile (radio sterilisation).

The implant is sealed in a film bag bearing a label for entering the implant data such as length and LOT number that can be put in the patient file. This film wrapping is the actual packaging.

The IMPLA Mini-Implant is intended to be used one time only!

If the packaging is damaged or not airtight, there is a risk of contamination with bacteria. The implant must not be used in this case.

General notes on implantation

Attention

All rotating instruments are subject to wear. Therefore, we cannot warrant the instruments' (pilot drills, extension drills, countersinks, cortical drills and thread cutters) functionality for more than 10 times of use, depending on the bone density. If the instruments are used more often, the manufacturer cannot guarantee a proper cutting efficiency and minimal temperature development any longer.

- Insufficient primary stability can be caused by poor preparation of the implant site or by insufficient bone quality. In this case, the loose implant should be removed and be replaced by an implant with a larger diameter. Implants should always be screwed into the bone all the way, even with their thread-free neck, when a covered healing process is considered (see above). Above all, a primary stability of the implant must be achieved.
- Implants should only be implanted when all prerequisites for a complication-free treatment and a longterm success are given. If there is any concern regarding the success of the implantation and/or the stability of the desired implant-based denture, a different solution without implantation should be favored.
- All rotating instruments can, in principle, be used several times. The user physician/expert alone decides under his/her own responsibility on the possibility of using the instrument again and the frequency of use on the basis of the respective case of use and the possible wear of the instruments. In case of doubt, the instruments should always be sorted out and replaced at an early stage.

Care of instruments and parts

- Used instruments/parts must always be disinfected before cleaning.
- Damaged instruments/parts must be disinfected, cleaned, and disposed of.

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Please consider our information on each material on page 30 – 37.

Disinfection

- Instruments/parts must be submerged in an adequate disinfection solution.
- Disinfection may be carried out manually or automatically. Please observe the manufacturer's instructions concerning dosing/concentration, residence time, temperature.
- In order to avoid damages, the instruments/parts may not have any contact during disinfection.
- It is highly recommended to wear protective gloves when working with contaminated instruments/parts (components).

Cleaning

- Only use adequate cleaning agents.
- Always observe the manufacturer's instructions concerning concentration of the agent, temperature, etc.
- It is highly recommended to wear protective gloves (DIN EN 455) when cleaning contaminated instruments/parts manually.
- In case of automatic cleaning and disinfection, the instruments/parts may not have any contact in order to avoid damages.
- Only instruments/parts made of the same material should be cleaned together.

- In case of heavy contamination of the instruments/parts, ultrasound cleaning techniques are recommended.
- After cleaning, rinse the instruments/parts thoroughly with water and dry them immediately.
- Damaged instruments/parts must be removed and disposed of separately.

Sterilization

- All surgical instruments and parts must be sterilized after cleaning.

Please consider our information on each material on page 30–37.

Storage

- The sterilized instruments/parts must be stored in a dry and dust-free atmosphere at room temperature.
- After the expiration of the maximum storage period, the instruments/parts must be sterilized again.

Dose 25 kGy/min

IMPLA locking screws (covering and healing caps) and gingiva formers

1. Indication

Surgical locking screws are used to immediately seal IMPLA implants after insertion and during submerged healing. The IMPLA gingiva formers can also be inserted during non-submerged healing or after successful submerged healing in order to form the gingiva surrounding the implant after the implant has healed.

2. Contraindications

None known

3. Notes

It is essential to read these instructions for use before using the IMPLA gingiva formers and locking screws.

When using the IMPLA gingiva formers and locking screws, the general rules for dental/surgical behaviour and indications, as well as the regulations for health and safety at work and the prevention of accidents must be taken into account. IMPLA gingiva formers and locking screws are only intended for medical/dental use in combination with the IMPLA implant system.

In the event of uncertainties in relation to the indication or type of application, use of the products should be ceased until all issues have been resolved. These instructions for IMPLA locking screws and gingivaformers are not sufficient to ensure proper use by handlers not experienced in implant-based prosthetics.

The user must be familiar with dental surgery and prosthetics, including diagnosis and preoperative planning and/or laboratory procedures. Since the use of this product is out of the control of the manufacturer, the latter shall not be liable in any way for damage caused. The responsibility rests exclusively with the user, or the attending physician.

When using IMPLA gingiva formers and locking screws, the following possible complications have been observed:

- The aspiration or swallowing of parts used in the patient's mouth
- Release or breakage of the screw thread due to incorrect torque (see recommended torque settings in Section 6 of the IMPLA instructions for use).

4. Precautionary measures

The following precautionary measures must always be taken before or during treatment.

- Before any procedures, ensure that all necessary parts, tools, and instruments are complete, functional, and available in the required amount.
- For your own safety, always wear suitable protective clothing.
- Position the patient so that the risk of aspirating the parts is minimised or excluded. Protect all parts used in the patient's mouth against aspiration and swallowing.
- Adhere to the specified torque for fixing prosthetic components (see recommended torque settings in Section 6 of the IMPLA instructions for use).

5. Side effects / Interactions

In very rare individual cases, allergies or sensitivities to the materials used to manufacture the prosthetic components (grade 5 titanium, etc.) cannot be ruled out.

In case of occlusal or approximal contact, the use of different alloy types in the same oral cavity may lead to galvanic irritation.

6. Use / insertion of locking screw / gingiva former

After fully inserting the implant in its final vertical position and loosening the insertion post, fit the surgical locking screw or gingiva former on the screwdriver (SW 1.2 mm), and screw into the implant by hand. The locking screws or gingiva formers are to be hand tightened (5–10 Ncm), taking into

account the primary stability of the implant. The gingiva former remains in for approx. 14 days to shape gingiva in place. The size of the gingiva former shall be selected according to the thickness of the gingiva, and the diameter and shape selected according to the prosthetic requirements.

Note: Gradually stretch to the desired diameter for very tight gingiva. Change to larger diameter every 5–7 days. Excessive pressure on the gingiva is to be avoided.

7. Cleaning / Sterilisation

The IMPLA locking screws (included in the implant packaging) and IMPLA gingiva formers are provided sterile in a suitable blister packaging.

8. Warning

The integrity of the product and sterile blister packaging must be checked before use. The product must not be used if the blister packaging is damaged.

IMPLA Abutments*

*The term „abutments“ is used here to represent ball connector, adhesive base, cone connector, SAE cone adapter, IMPLA Direct base.

Important notes / warnings:

IMPLA Abutments are a part of the IMPLA system and must only be used with the proper original IMPLA components. For the insertion of the abutments, only use the IMPLA instruments and tools. Otherwise, a damage of the original components cannot be ruled out.

Delivery:

All IMPLA Abutments are supplied in unsterile packaging. Therefore, the components must be sterilized accordingly before use.

Material:

IMPLA Abutments are manufactured from titanium grade 5 (Ti6Al4V ASTM 136).

IMPLA Scan abutments are made from PEEK (polyether ether keton).

Side-effects:

Allergies and sensitivities towards the alloy cannot be ruled out completely and may occur in very rare cases.

Interaction:

In case of occlusal contact, the presence of different alloy types in the patient's mouth may lead to galvanic reactions.

Safety advice:

Metal dust and ceramic dust are hazardous to the health. When working on these materials, always wear a protective mask and protective goggles. Use an extraction system with a fine particle air filter commonly used in dental practices.

Indications:

IMPLA ball abutments are used to anchor removable restorations to IMPLA implants. For this purpose, the prosthesis is ground out basally at a suitable area and a ball head matrix with an internal silicone O-ring is incorporated. The O-ring is subject to wear. It can be removed from the matrix and replaced, also for cleaning.

The IMPLA adhesive base is used to fix CAD constructed and CAM produced restorations on IMPLA implants. This requires a suitable design software as well as material blanks and a milling machine to produce the prosthetic structures, which can then be bonded to the adhesive base.

Cone connectors are intended for the production of cementable individual crowns and bridges.

SAE cone adapters are used for the production of bar abutments and bridges which are screw-retained in a tension-free method.

The IMPLA Direct base is intended for the construction of bar abutments or screw-retained bridges using the sprue technique, casting technique or laser technique.

Scan abutments are only intended for temporary use in the patient's mouth for the time of optical scanning using the intraoral scanner.

Contra-indications:

When selecting the patients, it is extremely important to ensure that they are able to maintain the necessary oral hygiene. In addition, the general rules and limitations in surgery and prosthetics for the selection and treatment of patients apply.

Procedure:

If the abutment does not bear a screw thread, carry out the final fixation on the intended implant with the vertical screw standard, art. no. 636649, with 20 Ncm, unless otherwise specified for the respective abutment type in separate instructions.

A few minutes after tightening the screw and at each recall (only for removable prosthetics) or at recalls with removal of the superstructure (for fixed prosthetics), it is advisable to check the seating of the screw again with the correct torque. If the structure is removed, a new vertical screw must be used when screwing it in again.

For the following cone connectors, the final fixation on the intended implant may only be performed with the vertical screw standard, article 636649 from batch 2020000804.

636198	Conical connector 3.3 mm / 15° – GH = 3 mm
636183	Conical connector 4.2 mm / 15° – GH = 3 mm
636184	Conical connector 5.3 mm / 15° – GH = 3 mm
636186	Conical connector 4.2 mm / 20° – GH = 3 mm
636187	Conical connector 5.3 mm / 20° – GH = 3 mm
636190	Conical connector 3.3 mm / 0° – GH = 1 mm
636191	Conical connector 4.2 mm / 0° – GH = 1 mm
636192	Conical connector 5.3 mm / 0° – GH = 1 mm
636193	Conical connector 3.3 mm / 0° – GH = 3 mm
636194	Conical connector 4.2 mm / 0° – GH = 3 mm
636195	Conical connector 5.3 mm / 0° – GH = 3 mm

For the following cone connectors, the final fixation on the intended implant may only be performed with the vertical screw short, article 638634.

636185	Conical connector 3.3 mm / 15° – GH = 1 mm
636196	Conical connector 4.2 mm / 15° – GH = 1 mm
636197	Conical connector 5.3 mm / 15° – GH = 1 mm
636188	Conical connector 3.3 mm / 20° – GH = 1 mm
636189	Conical connector 4.2 mm / 20° – GH = 1 mm

LOCATOR® ABUTMENT – Abutment System

The LOCATOR® Implant Attachment System includes: LOCATOR® Abutments, Retention Inserts, Denture Attachment Housing, Ancillary Processing Parts (including abutment analogs, processing spacer, block out spacer, parallelpost, impression coping, black processing male), and Tools (including LOCATOR® instrument, LOCATOR® adapter with connector for contra angle and corresponding adapter for torque ratchet).

DESCRIPTION: The LOCATOR® Implant Attachment System is a universal hinge, resilient attachment for endosseous implants in the mandible or maxilla in order to restore masticatory function. The attachment system allows for the prosthesis to be removed and replaced by the patient.

INDICATIONS FOR USE – The LOCATOR® Implant Attachment System is designed for use with overdentures or partial dentures, retained in whole or in part by endosseous implants in the mandible or maxilla.

CONTRAINDICATIONS – LOCATOR® Abutment: Not appropriate where a totally rigid connection is required. Use of a single implant with divergence of greater than 20 degrees from vertical is not recommended.

MRI SAFETY INFORMATION – A person with an implant or abutment device of the LOCATOR® Overdenture Implant System may be safely scanned under the following conditions. A scanning environment under Normal Operating Mode, with a Static Magnetic Field Strength (B_0) of 1.5 T or 3.0T, with a Maximum Spatial Field Gradient of 40 T/m (4,000 gauss/cm), and RF Excitation which is Circularly Polarized (CP). The Maximum Whole-Body specific absorption rate (SAR) is 2 W/kg and the Maximum Head SAR is 3.2 W/kg, both under Normal Operating Mode. The Scan Duration limit is 2 W/kg whole-body average SAR for 60 minutes of continuous radio frequency (RF), (a sequence or back-to-back series/scan without breaks). Failure to follow these conditions may result in injury. NOTE: The presence of this device may produce an image artifact that can extend as far as approximately 23.0 mm $\pm$ 0.5 mm from the device, as observed in a nonclinical setting. The $\pm$ 0.5 mm is the numerical value of the combined standard uncertainty.

STORAGE AND HANDLING – The LOCATOR® Implant Attachment System in its undamaged, original packaging is not subject to any special considerations for storage or handling (during transport and storage).

WARNINGS AND PRECAUTIONS – The abutments of LOCATOR® Implant Attachment System are made from the titanium alloy Ti-6Al-4V and have a Titanium Nitride (TiN) coating. Restorative components are made from Ti-6Al-4V, 303 Stainless Steel and Nylon. Patients should be evaluated for allergies or hypersensitivity to such materials prior to treatment. During the treatment process patients may be exposed to instruments, tools, and processing components which may contain various grades of stainless steel, polyethylene, nitrile, or silicone.

Product should be inspected for integrity prior to use. Product from damaged packaging should not be used on patients. In the event that the packaging is damaged, the damaged packaging with the product should be returned to the manufacturer and a replacement will be provided only if damage to packaging is caused by product shipment.

If the LOCATOR® Implant Abutment is subjected to inappropriate loading conditions, there may be a potential risk of metal fatigue.

As surgical instruments are susceptible to damage and wear, they should be inspected before each use. Markings should be visible and legible. Any reusable instrument should be replaced if damage or wear is present to ensure proper functionality. The number of uses will vary and depends on a variety of factors including but not limited to bone density encountered, handling, proper cleaning, autoclave exposure, and storage conditions (do not store tools or instruments wet). Over time, repeat sterilization may affect appearance and visibility of markings. When applicable to the surgical instrument, check the connection feature for wear to ensure the connection is not damaged.

Patient evaluation including the determination of the general health, oral hygiene habits and status, motivation towards good dental care, and anatomical acceptability prior to the placement of the implant attachments as part of restorative process, is critical. Thorough evaluation of the patient's medical

status and health history is mandatory. Treatment planning is vital to the success of the implant and prosthesis.

The use of each of these attachment systems require that the clinician be thoroughly familiar with the product and the method for its use and application. The clinician must also use reasonable judgment in deciding when and where to use the product.

SINGLE-USE DEVICES – The LOCATOR® Implant Attachment System components, with the exception of the tools and instruments, are single-use devices and are provided non-sterile. Single-use devices must not be reused or re-sterilized. Reuse of a single-use device may cause harm to the patient in the transfer of blood, tissue, or saliva that may contain infectious disease. Single use devices may not function as intended if re-sterilized and may result in an improper surgical procedure and lead to improper function or failure of the device.

LOCATOR® Males: The inadvertent re-use of LOCATOR® nylon males could cause loss of retention of the overdenture due to wear from previous use or damage during removal with the LOCATOR® Core Tool.

LOCATOR® Attachments: The inadvertent re-use of LOCATOR® Attachments could contain patient contamination build-up and subsequent wear of the retention inserts. This would result in improper fit and function which cause the loss of retention of the prosthesis.

MULTI-USE DEVICES – The surgical instruments and tools of the LOCATOR® Implant Attachment System are multi-use devices. Reusable tools and instruments must be cleaned and sterilized prior to reuse on patients

TOOLS: The LOCATOR® Tools are designed for multiple uses, and are provided NON-STERILE. Follow the instructions provided here within for proper sterilization of non-sterile components, and the instructions for cleaning and resterilization process of reusable components.

CLEANING AND STERILIZATION – The LOCATOR® Implant Attachment System, and other restorative components, instruments, and replacement LOCATOR® Attachments (sold separately), tools and instruments are supplied NON-STERILE, and should be sterilized prior to use on patients.

CLEANING – Reusable tools and instruments should be cleaned according to applicable instructions from the device manufacturer. **(i)** Disconnect or disassemble the instruments. **(ii)** Soak the instruments in enzymatic cleaning solution mixed according to manufacturer's instructions by completely submerging them for 20 minutes. Scrub the components using a soft-bristled, nylon brush until all soil has been removed. **(iii)** Remove the instruments from the enzymatic cleaning solution and rinse in tap water for a minimum of 3 minutes. Make sure to thoroughly flush internal holes/crevices of the instruments that have difficult to reach areas. **(iv)** Visually inspect instruments and tools for cleanliness and presence of residual debris. If additional cleaning is needed, place instruments in ultrasonic cleaner with enzymatic cleaning solution prepared according to manufacturer's instructions making sure that they are completely submerged, and sonicate for 10 minutes. **(v)** Remove the instruments from the ultrasonic cleaner, and rinse for 3 minutes making sure to thoroughly flush cleaning solution out of the holes/crevices and/or difficult to reach areas. **(vi)** Remove excess moisture from the instruments with a clean, absorbent, non-shedding wipe.

STERILIZATION – Tools and instruments provided non-sterile should be sterilized prior to use on patients. Reusable tools and instruments should be cleaned and sterilized prior to reuse on patients. The nylon males may be sterilized/disinfected using a liquid chemical sterilant as described below. Tools and instruments should be sterilized according to applicable instructions from the device manufacturer.

In order to ensure that the inserts are sterilized/disinfected, all microorganisms including *Clostridium sporogenes* and *Bacillus subtilis* spores are eliminated, the inserts must be soaked for a minimum of 3 hours in the liquid sterilant at room temperature. **Note:** An FDA approved liquid chemical sterilant for critical devices that are heat-sensitive and incompatible with sterilization methods such as steam and gas/vapor/plasma low temperature processes may be used following the manufacturer's directions for the sterilization of the device. Titanium abutments and stainless steel or other metal instruments may be sterilized by autoclave sterilization using the following parameters. For gravity cycle, place components in autoclave bag; and for pre-vacuum cycle, wrap the component with autoclave wrap

material and secure wrap with autoclave tape. Wrap the components using a wrap that is FDA-cleared for the indicated cycles.

Autoclave Sterilization Parameters are listed below:

Cycle Type	Description	Temperature	Exposure Time	Drying Time
Gravity	LOCATOR® Abutments, Tools and Instruments	132 °C / 270 °F	15 Min.	30 Min.
Pre-Vacuum		132 °C / 270 °F	4 Min.	20 Min.

Re-sterilizable instruments should be dried completely and stored in a clean and dry location at normal room temperature. Prior to instrument use, the exterior of any sterilized packaging should be inspected for integrity. Care must be exercised in the handling of wrapped or autoclave bagged instrument kits or instruments to prevent damage to the sterile barrier. If damage to the sterile barrier is observed, resterilization is recommended for reusable devices only. Single use devices should not be resterilized.

DISPOSAL – Dispose of used devices which pose a risk of infection according to facility clinical waste procedures and applicable local and state regulations.

PROSTHETIC PROCEDURES – Based on the results of the patient's pre-surgical assessment, the clinician should select and order the appropriate LOCATOR® Abutment based on the type of implant diameter, and tissue height.

It is imperative that all bone and soft tissue be removed from the superior aspect of the implant body to guarantee complete seating of the abutment.

Using a torque wrench, tighten the LOCATOR® Abutment to 30 Ncm. Use of higher torque values than recommended above could cause a fracture of the LOCATOR® Abutment.

Impression and Stone Model Fabrication

(i) With the LOCATOR® Abutments torqued in place, snap the impression copings on the abutments until they are seated firmly. (ii) Proceed by taking an impression. (iii) Remove the tray and snap an analog into each intaglio of the impression coping. (iv) Capture the abutment position in stone using standard methods for fabricating a laboratory stone model.

Prosthesis Fabrication

(i) Seat the LOCATOR® Retention Housings with the processing inserts on each of the abutments. (ii) Fabricate the prosthesis using standard laboratory techniques. (iii) When delivering the prosthesis, use the lowest retentive level insert to begin with and increase the retention if needed.

Denture Attachment Housing Pickup Technique (Optional) – (i) Place a block-out spacer around each abutment and press down.

(ii) Seat the LOCATOR® retention housings with the processing inserts on each of the abutments. (iii) Secure the retention housings to the prosthesis using auto-polymerizing or light cure acrylic or composite resin pickup technique.

Prosthesis Delivery

(i) Once the fit of the prosthesis is verified, remove the processing inserts from each retention housing using the core tool. (ii) Replace them with the lowest retention level inserts initially and increase the retention if needed. Firmly snap the prosthesis into place, ensuring that each insert is fully engaged onto each abutment.

HEALING PHASE – For delayed loading protocols: Relieve the denture to ensure the abutments are not in contact with any denture acrylic. A soft liner may be added to the denture to ensure patient comfort during the healing phase.

PATIENT CARE – Good oral hygiene is vital to attachment success. The patient should be made aware of the following:

(i) The LOCATOR® Abutments must be thoroughly cleaned each day to prevent plaque build-up and the patient should use a soft nylon bristle or end-tufted toothbrush with a non-abrasive toothpaste to clean the abutments. (ii) The coarse particles in abrasive toothpaste may scratch the surfaces of the abutments and cause plaque accumulation. (iii) An irrigation system is recommended to flush out debris from the inside of the LOCATOR® Inserts. (iv) The inserts are made of a soft plastic material (nylon) to allow the overdentures to be removed/replaced regularly. Plastic materials are subject to wear as part of normal use and may require replacement. (v) Bruxism wears the LOCATOR® attachments and may reduce the longevity of the retention inserts.

Patients should be instructed to maintain routine follow-up visits for hygiene and attachment function evaluation. Should a patient experience any discomfort or loss of retention of the overdenture, they should consult a dental professional.

Follow-up visits are recommended at 6 month intervals. Abutments must be re-tightened at follow-up visits to the torque specified above. Failure to re-tighten abutments could lead to screw loosening and abutment fracture. Patients should be examined for signs of inflammation around the implant abutments and for implant mobility.

Inserting and Removing the overdentures – The patient should be instructed on how to properly insert the overdenture. The patient should make sure they can feel that it is positioned over the abutments prior to applying pressure. The patient should use both hands and press down one each side and firmly snap the overdenture into place.

NOTE: The patient must not bite their overdentures into place as this force will result in improper wear of the abutments, including the retention inserts in the overdenture. Remove the overdenture by placing the thumbs under the edges of the overdenture flanges and pulling each side upward/downward simultaneously. Use of the tongue may aid in removal. Once removed, a thorough cleaning is recommended.

Cleaning your Implant Retained Denture – Instruct the patient to follow the protocol below to ensure the longevity of their overdenture. (i) Fill a washing basin with some warm water to prevent fracture of the overdenture. Apply non-abrasive toothpaste onto the soft bristle toothbrush and thoroughly clean every surface of the overdenture. (ii) Each night, remove overdenture and immerse in a cup of plain water.

Further Information – Traditional restorative protocols should be followed to process the attachments into the patient's overdenture. Standard overdenture care and maintenance should be followed in order to ensure the longevity of the restoration.

HSL Abutment Connectors

Description:

The IMPLA HSL Abutment Connector is used for the production of individualized, occlusally screw-fixed abutments in the cast-on technique. The abutment is composed of a base and a plastic tube (screw canal). The base is manufactured of a high-fusing cast-on alloy and the plastic tube burns out completely residue-free. During the cast-on procedure, the base is covered by the melt of a high-precious casting alloy or ceramic alloy. The metallic link is formed by diffusion of the components of both alloys.

Important notes / warnings:

IMPLA HSL Abutment Connectors are a part of the IMPLA system and must only be used with the proper original IMPLA components (IMPLA implants). For the insertion of the abutments, use only the original IMPLA instruments and tools. Otherwise, a damage of the original components cannot be ruled out.

Indications:

Single tooth gaps, partially or fully edentulous jaws, cement-fixed crowns or bridges, screw-fixed crowns.

Contra-indications:

Do not fix inappropriate abutments on IMPLA implants. Ensure an accurate selection of the abutment version (diameter of the connector and implant geometry). Manipulating the connector geometry will lead to an imprecise fit which will render the abutment useless.

Delivery:

All HSL Abutment Connectors are supplied in unsterile packaging. Therefore, the components must be sterilized by the practice staff accordingly before use.

Procedure:

Modelling	The plastic tube can be trimmed and modelled as customary in dental technology. Ensure the preservation of the connector geometry to the HSL-base and keep the connector geometry of the abutment absolutely clean to avoid cast-on in the area.
Investing	Ensure that the investment remains bubble-free. Only use phosphate-bonded investment materials when casting high-precious metals.
Preheating	In order to allow for a complete cast-on, please adhere to the instructions of the preheating furnace as well as to the preheating temperatures and holding times prescribed by the manufacturer of the alloy.
Casting	Please adhere to the instructions of the casting machine. To prevent the IMPLA HSL Abutment Connectors from melting, please ensure that the casting temperature of the ceramic alloy you use is significantly lower than 1400°C (solidus temperature).
Cooling	We recommend letting the casting cylinder cool to room temperature to avoid tension in the cast object.
Divesting	To avoid a deformation of the object, especially of the connector geometry to the implant, only blast with beads at low pressure. Afterwards, the object can be stripped in an appropriate stripping solution.
Finishing	We recommend using ceramic grinders and hard-metal mills. Adhere to industrial safety regulations (protective goggles, extraction system, respiratory protection)!
Veneering	Please note the coefficient of thermal expansion of each veneering material to prevent cracking.

The abutments are fixed with the corresponding Aesthura vertical screw at 20 Ncm.

A few minutes after tightening the screw and at each recall (only for removable prosthetics) or at recalls with removal of the superstructure (for fixed prosthetics), it is advisable to check the seating of the screw again with the correct torque. If the structure is removed, a new vertical screw must be used when screwing it in again.

Composition of the IMPLA HSL-Base:

Gold	60 %
Platinum	24.9 %
Palladium	15 %
Iridium	< 1 %

Technical data:

Density	18.1 g/cm ³
Melting interval	1350 – 1460°C
Vickers hardness HV5/30	150 – 205 – 230

Side effects:

Allergies and sensitivities towards the alloy cannot be ruled out completely and may occur in very rare cases.


Interaction:

In case of occlusal contact, the presence of different alloy types in the patient's mouth may lead to galvanic reactions.

Multi-Unit Abutments

Description:

The IMPLA Multi-Unit Abutment Concept facilitates improved usage of the available bone structure. By means of an angulated insertion of the terminal implants, the available bone is used in distal direction which improves the stability. With the according angulation, the Multi-Unit Abutments ensure a uniform insertion direction of the prosthetics to facilitate a tension-free screw fixation. For this purpose, the IMPLA Multi-Unit Abutment Concept contains straight and angulated components.

Components Multi-Unit System

Abutments with the diameters 3.3 mm, 4.2 mm, 5.3 mm, each with two gingiva heights of 1 mm or 1.5 mm and 3 mm, and angulations of 0°, 20° and 30° (not 5.3 mm).

Accessories IMPLA Multi-Unit System

Impression system with impression aid and fixing screw as well as the according lab implant.

Gingiva sleeve for conditioning of the mucosa as well as the appropriate metal or plastic sleeve for the prosthetics with suitable screw for secondary screw fixation.

Important notes / warnings:

IMPLA Multi-Unit Abutments are a part of the IMPLA system and must only be used with the proper original IMPLA components (IMPLA implants). For the insertion of the abutments, use only the IMPLA instruments and tools. Otherwise, a damage of the original components cannot be ruled out.

Indications:

The IMPLA Multi-Unit Abutments are exclusively meant for a use with occlusal screw-fixed prosthetics (e. g. bridges, hybrid prostheses and cast bars).

A minimal implant length of 11.5 mm in the distal position is indicated for the use of IMPLA Multi-Unit Abutments.

The minimal number of implants for a complete reconstruction is four in the mandible and six in the maxilla.

A sufficient bone quality is a basic requirement for implantation, as it ensures the primary stability. The proper torque is between 30 and 40 Ncm.

The distal implants must be inserted with an angulation that corresponds to the angulation of the abutments inserted later, i.e. **0°, 20° or 30°**.

Presurgical planning as well as a comprehensive medical history and diagnosis (including radiological examination) are absolutely necessary.

In principle, the general rules and instructions for surgical and prosthetical procedures in dental implantology apply.

Contra-indications:

When selecting the patients, it is extremely important to ensure that they are able to maintain the necessary oral hygiene. In addition, the general rules and limitations in surgery and prosthetics for the selection and treatment of patients apply.

Delivery:

All components of the IMPLA Multi-Unit Abutment Concept are supplied in unsterile packaging. Therefore, the abutments must be sterilized accordingly before use.

Procedure:

Based upon presurgical planning and following the drill protocol for IMPLA implants, a pilot drilling for a mesial implant without angulation is performed first.

Subsequently, all other implants can be aligned to this first drilling with help of a drilling template. The drilling template improves the spacial orientation and presets the angulation of the distal implants. The drill is guided through a sleeve in the drilling template into the bone in the exact angle intended.

After the cavities have been prepared in accordance with the IMPLA drilling protocol, the implants can be inserted.

The prosthetic restoration can be carried out either in two stages or by immediate load.

In case of a two-stage restoration, we recommend to cover the implant with a surgical cover screw (covered healing) or with a gingiva former for open healing.

A primary stability of 30–40 Ncm is absolutely necessary for the immediate load approach.

Impressions are taken with help of the according impression components of the IMPLA Multi-Unit Abutment Concept.

Depending on the previously selected implant angulation (**0°, 20° or 30°**), the mesially positioned implants are then restored with a straight abutment and the distally positioned implants with an angulated abutment. The abutments are fixed with the appropriate vertical screw with a torque of 20 Ncm.

A few minutes after tightening the screw and at each recall (only for removable prosthetics) or at recalls with removal of the superstructure (for fixed prosthetics), it is advisable to check the seating of the screw again with the correct torque. If the structure is removed, a new vertical screw must be used when screwing it in again.

Finally, the prosthetic construction is set upon the abutments tension-free and is fixed with a prosthetic screw with a torque of 15 Ncm.

Material:

The IMPLA Multi-Unit Abutments are made of titanium grade 5 (Ti6Al4V ASTM 136).

Side effects:

Allergies and sensitivities towards the alloy cannot be ruled out completely and may occur in very rare cases.

Interaction:

In case of occlusal contact, the presence of different alloy types in the patient's mouth may lead to galvanic reactions.

IMPLA PreFace Abutments

Description:

The IMPLA PreFace Abutment is finished with an appropriate milling system with help of the original PreFace Abutment holder. The PreFace Abutment holder serves to fix the PreFace Abutment inside the milling system in the appropriate working position. The PreFace abutments are held in the PreFace abutment holder on the front side, making sure that they are seated exactly in the end position. The PreFace abutment holders are available for a wide range of milling machines.

Important notes / warnings:

IMPLA PreFace abutments are part of the IMPLA system and must only be used with original components of the system (IMPLA implants). Process the PreFace abutments only with the corresponding suitable abutment holders.

Indications:

The IMPLA PreFace Abutments are used for the production of individual titanium abutments on IMPLA implants in CAD/CAM milling machines. The individual abutments can be combined with crowns or bridges to reconstruct function and aesthetics.

Contraindications:

The IMPLA PreFace Abutments must only be combined with IMPLA implants. Manipulating the connector geometry will lead to an imprecise fit which will render the abutment useless.

Delivery:

All components of the IMPLA PreFace Abutment concept are supplied in unsterile packaging. Therefore, the components must be sterilized accordingly before use by the practice staff.

Procedure:

The shape of the PreFace Abutment can be adjusted individually to the present anatomic conditions. The production of the PreFace Abutment is done digitally with an appropriate CAD software. The basic rules must be adhered to.

Design the emergence profile gently. Any further expansion of the gingiva than the expansion obtained with the gingiva former should always be discussed with the attending dentist beforehand. The gingiva can be expanded gradually by individually designed temporary abutments.

In general: The circular notch should be positioned slightly towards the subgingival direction in the vestibular area and orally on the gingiva line. This way, the cement between coping and abutment is easy to remove and the transition is not visible from vestibular direction. The minimal wall thickness of the abutments should not fall below 0.5 mm (screw hole to the outer surface of the abutment). The diameter of the emergence profile should not exceed 13 mm, while the height of the emergence profile should not exceed 6 mm. The total length of the PreFace Abutment should not exceed 15 mm, the maximum angle of the PreFace Abutment should not exceed 30°. The taper of the PreFace Abutment should not exceed 6°. The PreFace Abutment should be rounded occlusally and should taper off. Avoid sharp edges.

The recommended torque for the fixation of IMPLA PreFace Abutments is 20 Ncm.

A few minutes after tightening the screw and at each recall (only for removable prosthetics) or at recalls with removal of the superstructure (for fixed prosthetics), it is advisable to check the seating of the screw again with the correct torque. If the structure is removed, a new vertical screw must be used when screwing it in again.

Finishing:

In case that further work steps are necessary after shaping the PreFace Abutment in the CAD/CAM milling machine, please note the following:

To mill PreFace Abutments, use appropriate, flawless mills with only little pressure. The minimal wall thickness must not fall below 0.5 mm. Avoid ridges and edges. We recommend to fix the PreFace Abutments in the lab analog to protect the connection geometry.

Scan body:

A scan abutment is used to determine the exact position of the implant during the scanning procedure. The scan abutment's diameter must correspond to the diameter of the selected implant.

* PreFace is a registered trademark of Medentika GmbH

IMPLA Aesthura Abutments

Description:

The IMPLA Aesthura abutment is used to fabricate occlusally screw-retained restorations. It is screwed onto the intended implants with the vertical screw, art. no. 638665.

Important information / warning:

The Aesthura Abutment must not be altered by measures of any kind.

IMPLA Aesthura Abutments are part of the IMPLA system and must only be used with original components of the system (IMPLA implants). For the insertion of the abutments, use only the IMPLA instruments and tools. Otherwise, a damage of the original components cannot be ruled out.

Indications:

Individual gaps, partially or completely edentulous jaws, screw-fixed individual tooth restorations.

Contra-indications:

When selecting the patients, it is extremely important to ensure that they are able to maintain the necessary oral hygiene. In addition, the general rules and limitations in surgery and prosthetics for the selection and treatment of patients apply.

Delivery:

All IMPLA Abutments are supplied in unsterile packaging. Therefore, the components must be sterilized accordingly before use by the practice staff.

Material:

The IMPLA Aesthura Abutments are made of grade 5 titanium (Ti6Al4V ASTM 136).

Side effects:

Allergies and sensitivities towards the alloy cannot be ruled out completely and may occur in very rare cases.

Interaction:

In case of occlusal contact, the presence of different alloy types in the patient's mouth may lead to galvanic reactions.

Procedure:

The abutments are fixed with the corresponding Aesthura vertical screw at 20 Ncm.

A few minutes after tightening the screw and at each recall (only for removable prosthetics) or at recalls with removal of the superstructure (for fixed prosthetics), it is advisable to check the seating of the screw again with the correct torque. If the structure is removed, a new vertical screw must be used when screwing it in again.

Disinfection, cleaning and sterilization of surgical and prosthetic tools and instruments of the IMPLA system

Following surgical and prosthetic instruments/tools of the IMPLA system are composed of surgical steel/ titanium/PEEK/PP and are to be treated accordingly:

X-ray spheres, drills, drill extensions, countersinks, bone thread cutters, ratchet, screw drivers, IPK, depth stop, gingiva cutter, contra angle, flat wrench, sleeve holder, tray, carrier and guide sleeves, depth gauge, tweezers, sealing and healing cap, gingiva former, abutments, scan abutments, impression aid, horizontal rosette, vertical screw, bar screw, bone screws

Please observe!	Neutralize and dry the instruments thoroughly in order to avoid precipitation on the instruments to be sterilized. It is recommended to autoclave unmixed materials only.
Disinfection	Always wear protective gloves when disinfecting components suited to this treatment. Avoid drying of sticking residues of blood and/or tissue. Always use an appropriate disinfection bath and dose according to the manufacturer's instruction.
Cleaning	The components can be cleaned with all cleaning agents suitable for reprocessing surgical instruments and are resistant to any cleaning process, e.g. ultrasonic bath. Please avoid the use of metal brushes and abrasives. After cleaning, subsequently clean pH-neutrally and rinse off the cleaner and particle residues with demineralized water. Dry immediately.
Sterilization	The components are resistant to any kind of sterilization process. Make sure that all deposits have been removed by cleaning. The parts may not have any contact during sterilization. The temperature should not exceed 176°C (248°F). Please observe the instructions of the autoclave manufacturer.

Instructions for the practice:

1. Do not let surgical residues (blood, tissue, etc.) dry on the instruments/parts. Submerge the instruments/parts in disinfective agent directly after the surgery.
2. Clean the instruments/parts only with a nylon brush. (Do not use metal brushes, such as brass brushes, etc.)
3. Rinse off disinfectant residue thoroughly with clear water.
4. Dry instruments/parts directly after cleaning and rinsing.
5. Never sterilize instruments made of different materials (e.g. ratchet) in the same cycle!
6. Please observe the instructions for use of the cleaning agents, as suitability for all instruments is limited (see material groups).

1 Use of the torque ratchet:

The torque ratchet is build up of the ratchet head (Fig. 1-b), the ratchet handle (Fig. 1-c) and the setting screw (Fig. 1-e). Inside the ratchet head, there is a tool fitting (Fig. 1-a).


Fig. 1

The ratchet head is lockable to lock the torque limit (see 1.1 below). If the ratchet head is not locked, the torque can be adjusted infinitely between 10 and 30 Ncm with the setting screw (see 1.2).

Before its first use and after each use in a patient case, the ratchet must be disassembled, cleaned and dried. **It must be sterilized directly before each use** (see section 3).


Fig. 2


Fig. 3

Fig. 4


Fig. 5

1.1 Adjusting the ratchet

In the factory setting, the ratchet head is not locked. If you would like to lock it, turn the setting screw to less than 10 Ncm (Fig. 4), to relieve all pressure on the spring.

Subsequently, pull the handle from the ratchet head (Fig. 2). Then, turn the handle clockwise by 90° and let it lock into place in the ratchet head. Finally, preset a torque of 20 Ncm to ensure that the ratchet remains locked.

Attention: In the fixed condition, the maximum load is 100 Ncm!

1.2 Setting the torque function

First, ensure that the ratchet is not locked. Turn the setting screw until you have reached the necessary torque. The correct torque is set once the according marking ring at the end of the handle is just still visible (Fig. 4). The ratchet head will now bend once the preset torque has been reached (Fig. 5). Afterwards, do not keep on turning, as this could lead to exceeding the preset torque significantly.

Important notice:

Whenever the ratchet is not in use, turn back the setting screw to relieve all pressure from the spring.


1.3 Change tool fitting

(See 2.1 [1.] and 2.2 [5.])

2 Maintenance of the torque ratchet

2.1 Disassembly

1. The tool fitting (Fig. 6-a) is removed from the ratchet head (Fig. 6-b) by pulling on the two pin ends (arrows) towards the ratchet handle. Please note that the tool fitting can only be pulled from the ratchet head towards one side. **Do not use force.**
 2. Completely screw off the setting screw (Fig. 6-e) from the ratchet handle (Fig. 6-c).
 3. Pull off the ratchet handle (Fig. 6-c) from the ratchet head (Fig. 6-b) and remove the spring.
 4. Disinfect, clean and dry all parts (see section 3).
- Afterwards and **before the sterilization**, the ratchet is reassembled.

2.2 Reassembly

1. Before reassembling the ratchet, all parts must be disinfected, cleaned and dried completely. Cover the areas marked with ● with oil (see 3.4).
2. Push the ratchet handle (Fig. 6-c) over the threaded bolt into the ratchet head (Fig. 6-b).
3. Push the spring (Fig. 6-d) into the ratchet handle (Fig. 6-c).
4. Screw the setting screw into the ratchet handle (Fig. 6-e) until the resistance of the spring becomes evident.
5. The tool fitting (Fig. 6-a) is set back into the ratchet head (Fig. 6-b) by pulling the two pin ends (arrows) towards the ratchet handle. Please consider that the tool fitting can only be set into the ratchet head from one side. **Do not use force.**

Important information:

Always relieve pressure on the spring when ratchet is not in use. Screw off the setting screw until the pressure on the spring is relieved completely.

3 Reconditioning

3.1 Disinfection

Before first use and immediately after each use, all ratchet parts must be disinfected in a suitable disinfectant. Adhere to the instructions of the disinfectant's manufacturer.

Chemical as well as thermal disinfection techniques are suitable, if they come up to the requirements of medical instrument disinfection.

3.2 Cleaning

Immediately after disinfection, all parts of the ratchet must be cleaned with a suitable cleaning agent. Adhere to the instructions of the cleaning agent's manufacturer. Otherwise, impurities will dry off onto the instrument and make cleaning difficult.

In order to avoid corrosion and discoloration, closely adhere to the instructions on concentration, submergence time and temperature given by the manufacturer of the cleaning agent.

Remove any deposits, impurities, etc. thoroughly with a soft brush (do not use metal brushes, scouring agents, wire wool or too much pressure).

Please note:

When disinfecting or cleaning the device with an ultrasound unit, please adhere to the instructions of the unit's manufacturer.

3.3 Rinsing and drying

Directly after disinfection, cleaning or ultrasound treatment, the individual parts must be rinsed carefully with clear desalinated or distilled water (do not use tap water). Afterwards, dry all parts thoroughly to avoid water stains.

If the individual parts are not rinsed, residue is not removed and/or the parts are not dried completely, this could lead to discoloration or corrosion of the device during sterilization.

3.4 Care and checking

Parts that have joints must be treated with a sterilizable oil spray (paraffin-based) after all cleaning and disinfection procedures. Such a spray will reduce metal-on-metal-friction, avoids fretting corrosion and reduces the danger of oxidation during heavy use.

3.5 Sterilization

Only sterilize carefully cleaned, dried devices without any type of corrosion.

Badly cleaned and insufficiently dried parts bear the danger of burn-in or corrosion damage.

If corroded devices are sterilized with immaculate instruments, flash rust can cause damages on the immaculate instruments.

Please note the instructions of the manufacturer.

Note:

Only sterilize the torque ratchet directly before use.


3.6 Storage

Important information regarding the storage of the ratchet:

The spring of the ratchet should be fully relieved immediately after use to ensure proper use in the long term. Only store the ratchet in a relieved state.

Cleaning and sterilization of the ratchet

Screw off turn-knob, remove spring and handle.

To clean, pull back the interlock to unlock the receptacle

*Sterilization of the ratched only in disassembled state, without the spring.
Attention: Autoclave spring separately!*

Application and safety instructions

All rotating instruments and parts for medical use (dental implantology, surgery) were developed for their specific application. For that reason, incorrect use may cause premature wear of the instruments and endanger the patients.

Rotating instruments for the medical field must only be used by physicians and other experts who are intensively familiar with the use of these instruments due to their training and experience.

Correct application

- Take care to use only technically and hygienically flawless, well-maintained and cleaned contra-angle handpieces.
- The instruments must be inserted correctly.
- The instrument must revolve at operating speed before being placed on the object (see pages 4, 8 and 12).
- Avoid wedging or levering the instrument as this may increase the risk of fracture.
- Thermal damage through the use of rotating instruments must be avoided (we recommend the slow-speed mode with sufficient cooling). Please see Surgical procedure or Procedure on pages 4, 5, 8 and 12 of this manual.

Pressing forces

- Excessive pressing forces must be avoided at all times. Excessive pressing forces may damage the working sections of the instruments or fracture their blades. This also results in an increased generation of heat.

Cooling

- To prevent excessive heat generation during preparation, provide for adequate air or water cooling (at least 50 ml/min).
- Inadequate cooling may lead to irreversible damages of the bone and/or tissue.

Manufacturer's information on rotating surgical instruments

Product:

The present manufacturer's information applies to all instruments that are used for surgical application. These are rotating instruments made of stainless steel.

Instruments delivered in non-sterile condition must be prepared prior to first use.

Limited number of reprocessing cycles:

The end of product life is basically determined by wear and damage from use. If necessary, known restrictions on the frequency of use of instruments must be observed (see pages 5, 8 and 13). Frequent reprocessing does not affect the performance of these instruments.

Work station:

Hygienic precautions according to the provisions valid in your country.

Storage and transport:

Place instruments in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). It is recommended to process the instruments within one hour of use at the very latest. The instruments should be transported to the reprocessing site in the disinfection tank. Endodontic instruments can be transported in a special interim support equipped with a foam insert drenched in disinfecting solution.

Standardized manual reprocessing for rotating instruments

Equipment used:

- Nylon brush
- Suitable detergent/disinfectant for rotary instruments with proven disinfecting effect
- Ultrasonic bath

Reprocessing:

1. Rinse off surface contamination under running water. Remove stubborn contamination with a nylon brush under running water, turning the instrument constantly.
2. Rinse instrument thoroughly with running water.
3. Place the instruments in a suitable sieve into the ultrasonic device filled with detergent/disinfectant.
4. Observe the instructions of the manufacturer regarding concentration and immersion time for chemical cleaning/disinfection in the ultrasonic device. Be sure to observe the full correct immersion time which does not start until the last instrument has been placed into the ultrasonic device. Attention: Do not exceed 45°C (113°F) (risk of protein coagulation)!
5. On completion of the immersion time, rinse instruments thoroughly with suitable water (preferably with demineralised water to avoid residues).
6. Dry instruments (according to RKI recommendations, preferably with compressed air).
7. Visual examination with a suitable magnifying device to ensure that the instrument is clean and undamaged (experience has shown that a magnification factor of 8 permits a visual examination). If there are still visible residues of contamination, repeat the cleaning and chemical disinfecting process until no visible contamination is left.

Control and functional test:

Instruments showing the following defects are to be discarded immediately:

- Blunt and chipped blades
- Deformations (e. g. bent/twisted/fractured instruments)
- Corroded surfaces

Packing:

Make sure that the packaging is suitable for the instrument and the chosen method of sterilization. Single pack: The packaging must be large enough to ensure that there is no pressure on the seal. In the set: Place instruments onto the tray provided or onto universal sterilization trays. The instruments must be protected. Use an appropriate method to pack the tray. Instruments with limited times of use are to be marked accordingly.

Sterilization:

Steam sterilization using a vacuum process at 134°C (273°F) in a device that complies with the provisions of DIN EN 13060; with validated processes.

- Fractionated pre-vacuum (type B) respectively simplified pre-vacuum (type S).
- Sterilization temperature: 134°C (273°F)
- Hold time: 5 minutes minimum (full cycle)
- Drying time: 10 minutes minimum

In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.

Transport and storage:

The packed sterile goods must be protected from dust, moisture and recontamination during transport and storage.

General note:

Observe the legal provisions regarding the reprocessing of medical products valid in your country (e. g. https://www.rki.de/EN/Home/homepage_node.html).

The manufacturer confirms that the above detailed reprocessing methods are **suitable** for preparing the above named instrument group to enable their reuse. The person carrying out the reprocessing is responsible for ensuring that the applied method is carried out with appropriate equipment, materials and trained personnel at the reprocessing site and that it actually achieves the desired result. To guarantee this, routine controls of the validated mechanical and/or manual preparation methods are necessary. Any deviation from the above detailed process (e. g. use of different chemicals) must be carefully checked by the operator to ensure effectiveness and to avoid possible adverse consequences.

Information on storage, disinfection, cleaning procedures and sterilization of rotating instruments

- Rotating instruments must be disinfected, cleaned and – whenever the instrument material permits it and it seems suitable – sterilized prior to first use on patients and immediately after use.
- Used instruments must always be disinfected before cleaning.
- Before first use, the rotating instrument should be kept in its original packaging and be stored at room temperature, protected against dust and moisture.
- All instruments and parts are packaged non-sterile and must be sterilized before use according to their application. Cleaning and sterilisation not in accordance with the instruction may lead to microbiological contamination and infection of the patient.
- All materials can be sterilized in an autoclave according to DIN ISO 17665 „Sterilization of health care products - Moist heat“ (e. g. at 134°C [273°F] with a minimum holding time of 3 min. or 121°C [249°F]/15 min.). Please be sure to follow the manufacturer's operating instructions when sterilizing.
- Before sterilization, the instruments must be cleaned and disinfected properly.
- Disinfection may be carried out manually or by machine.
- In case of severe contamination of the instrument, ultrasound cleaning techniques are recommended.
- After cleaning rinse the instruments thoroughly with water and dry them immediately.
- Contact between the instruments during cleaning and disinfection procedures may cause damage and should be avoided.
- It is highly recommended to wear protective gloves when working with contaminated instruments.
- Thermal disinfection procedures are generally not suitable for rotating instruments.
- Avoid sterilization temperatures exceeding 180°C (356°F). Exceeding this temperature reduces the hardness of the working part and reduces the instrument's service life.
- All sterilized instrument containers should be labeled with an indicator strip bearing the date of the sterilization as well as the expiry date.
- Non-rustproof instruments may corrode. This may cause discolorations and reduce the service life. Therefore, it is recommended to use anti-corrosion disinfectants and cleaning agents in these cases.
- The instructions for use, reaction time, and suitability of disinfectants and cleaning agents for certain types of instruments are provided by the manufacturer's instructions. They must be observed.
- Contact of the instruments with hydrogen peroxide (H₂O₂) must be avoided.
- The instruments should be stored in appropriate, hygienically well-kept containers. This also applies to sterilized instruments. When stored, the instruments should be protected against dust, humidity, and recontamination. Observe the maximum storage period.
- Broken-out cutting edges cause vibrations and high contact forces and thus lead to broken-out preparation edges and rough surfaces. They should be discarded immediately.
- Bent or non-concentric instruments should be discarded immediately.