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Sebond pink

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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Sebond pink

Intended use:

Framework covering materials are light-curing or chemically curing composite-based materials for covering different framework materials in preparation for subsequent veneering with veneering composite or denture resin.

Product description

Sebond pink is a gingiva colored, light-curing paste opaque with integrated bonding component designed for covering framework retentions.

General product description:

- Light-curing one-component material suited for color coating framework materials (conventional dental metals and alloys). The material is light-curing and is to be polymerized in a light-curing unit. Sebond pink contains an additional bonding component which assures its bond to the framework material. The use of a supplementary bonding agent is therefore unnecessary.
- The material is applied to the retentions of the frameworks to cover them in the gingival shade before they are completed with acrylic.
- Sebond pink is available in the shades "light" and "dark" .

Composition

- Silicon dioxide, glass powder, monofunctional aliphatic urethane acrylate, aliphatic urethane dimethacrylate, pigments, initiators
- Filler content: 45 % by weight of inorganic filler (0.005 – 3.0 µm)

Indications

Covering framework retentions (conventional dental metals and alloys) when adding a gingiva restoration.

Contraindications

- Do not use the product in case of known allergies against one of the components.

Interactions with other agents

Phenolic substances (e. g. Eugenol) inhibit polymerization. There-fore, do not use materials which contain such substances.

Application/processing

Processing time of the masses: 1 – 3 min., depending on lighting conditions.

1. Preparatory work

The framework is fabricated according to the manufacturer's specifications. The surface is sandblasted with aluminum oxide 50 – 110 µm, 2 bar and blown off with oil-free compressed air.

2. Procedure

Sebond pink is applied in several layers (min. 2 layers), wash-like. Each layer is polymerized separately (see table Polymerization times).

3. Subsequent procedures

The acrylic can be cast on subsequently.

Troubleshooting / FAQ list

- The framework surfaces must be cleaned of grease layers and polishing residues.
- Sebond pink must not be mixed with opaque liquids or opaque powders
- Do not use intermediate polymerization equipment for poly-merization.

Problem	Cause	Remedy
does not cure	applied layer is too thick	max. layer thickness 2 mm
	insufficient polymerization	observe poly-merization times – check lamp, replace lamp if necessary – extending the polymerization times does not result in greater depth of cure – do not use an intermediate poly-merization lamp
chipping in the opaque	opaquer not poly-merized sufficiently/ applied layer too thick	only apply wash-like respectively in very thin layers incorrect polymer-ization times/ check lamp/ replace lamp if necessary
material too firm inside the syringe	syringe stored at a temperature below 10 °C (50 °F)	store material at room temperature or let it reach room temperature before use
surface greasy	insufficient polymerization	observe poly-merization times
		check unit/ maintain regularly

General information

A suitable light polymerization unit with an emission spectrum in the range of at least 310 – 500 nm must be used. The required physical properties can only be achieved with perfect lamps. There-fore, a regular check of the light intensity according to the manu-facturer's specifications is required. Do not use intermediate poly-merization units.

Storage

- Storage temperature 10 – 25 °C (50 – 77 °F).
- Close syringe tightly immediately after use

Shelf life

The maximum shelf life is printed on the label of each syringe. Do not use after expiration date.

Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-indications / interactions

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

Note

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Polymerisation times for Sebond pink:

Light curing device	Time
Spekra LED	1 min.
Spekra 2000	3 min.
HiLite/ Uni-XS	90 sec.
Labolight LV-II/ III	1 min.
Solidilite	1 min.

CE 0297

SCHÜTZ DENTAL
Micerium Group

Schütz Dental GmbH
Dieselstr. 5-6 · 61191 Rosbach/Germany
Tel.: + 49 (0) 6003 814-0 · Fax: + 49 (0) 6003 814-906
info@schuetz-dental.de · www.schuetz-dental.de