dialog Vario

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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Instructions for Use

These instructions for use comprise the dentin masses, incisal masses and effect masses of the products dialog Vario and dialog Vario Flow.

Purpose

Veneering resins are light-curing composites of different consistencies for coating (veneering) various framework materials for the design of dental restorations such as e.g. crowns and bridges.

Product description

dialog Vario and dialog Vario Flow are used for veneering various framework materials. In this way, an aesthetic, tooth-shaped appearance of the restoration is achieved.

Different consistencies allow users to use different modeling techniques for the design. In addition, a wide range of shades facilitates a perfect imitation of shape and color.

Users

For use by dental technicians in a dental laboratory.

Composition dialog Vario

- Glass powder, splinter polymer, urethane dimethacrylate, tetramethylene dimethacrylate, Bis-GMA, silicon dioxide, pigments, initiators
- filler content: 67% (by weight) anorganic filler (0.005 - 3.0 µm)

dialog Vario Occlusal

- Glass powder, urethane dimethacrylate, silicon dioxide, Bis-GMA, tetramethylene dimethacrylate, pigments, initiators
- filler content: 75% (by weight) anorganic filler (0.005 - 3.0 µm)

dialog Vario Flow

- Glass powder, urethane dimethacrylate, tetramethylene dimethacrylate, silicon dioxide, pigments, initiators
- filler content: 56% (by weight) anorganic filler (0.005 - 3.0 µm)

Indications

Fixed restorations

With framework support

- Veneering crowns, bridges, adhesive bonding bridges, implant constructions, supra constructions as well as long-term temporaries

Without framework support

- Production of inlays, onlays, indirect veneers and framework-free full veneers of individual tooth restorations

Removable restorations

With framework support

- Veneering conic and telescopic crowns as well as attachment exterior parts

- Individualization of prefabricated acrylic teeth
- Repairs, cosmetic corrections and characterization of existing restorations

Contra-indications

Do not use in case of a known allergy against one of the components.

Interaction with other agents

- Phenolic substances (e.g. eugenol) inhibit polymerization. Do not use any materials which contain such substances.
- Generally known cross-reactions or interactions of the medical device with other materials already in the mouth must be taken into account by the dentist when using the product.

Framework design and preparation

- The frameworks are modeled, cast or CAD/CAM milled and finished as usual in accordance with the applicable dental technology guidelines.
- Conventional dental metals and dental alloys, PMMA, PEEK and zirconium dioxide can be used as framework materials. All materials must be processed and prepared for further processing according to the manufacturer's instructions.
- The frameworks must be conditioned before the composites are used. Depending to the framework material used, this framework conditioning is carried out with a suitable bonding agent in accordance with the respective manufacturer's instructions.
- The framework surfaces can also be covered with a paste opaque. Please refer to the manufacturer's instructions.

Before final polymerization, the entire veneering surface can be covered with experts gel. This prevents the formation of a new dispersion layer and facilitates finishing.

Veneering of fixed or removable restorations with framework support

(i. e. veneering of crowns, bridges, adhesive bonding bridges, implant constructions, supra constructions, long-term temporaries, conic crowns, telescopic crowns, attachment exterior parts, as well as individualization of prefabricated acrylic teeth, repairs, cosmetic corrections and characterization of existing restorations)

- **Note:** The incisal and effect masses must **not** be applied **directly** to a framework, a bonder surface or an opaque surface. In a first step, a layer of dentin material is therefore applied to the framework after the bonding agent and/or opaquer has been applied. This serves as an elastic or flexible intermediate layer and absorbs forces and stresses. After layering the incisal materials or effect materials, the final polymerization is carried out, see polymerization table.
- After applying the according bonder and opaquer, remove the desired **dentin material**, paste or flow from the syringe, apply to the area to be veneered and shape into the desired form using a brush or spatula.
- The maximum layer thickness is 2 mm. If the area to be completed exceeds this thickness, intermediate polymerization must be carried

Color combination table

Color according to Vita®**	Dentin	Incisal
A1	A1/B1	S 58
A2	AZ/B2	S 58
A3	A3	S 59
A3,5	A3,5	S 59
A4	A4	S 60
B1	A1/B1	S 57
B2	AZ/B2	S 58
B3	B3/B4	S 59
B4	B3/B4	S 59
C1	C1/D2/D3/D4	S 59
C2	C2/C3/C4	S 59
C3	C2/C3/C4	S 59

C4	C2/C3/C4	S 60
D2	C1/D2/D3/D4	S 60
D3	C1/D2/D3/D4	S 59
D4	C1/D2/D3/D4	S 60

* Vita is a registered trademark of Vita Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen/Germany.

Recommendation for layering

B2: On the edges, draw incisal all the way down to the shoulder area, but cover the shoulder area only with a relatively thin layer of incisal material.

C1: Completely cover tooth body with incisal material.

C4: Modify dentin C3 with khaki/blue mixture. On the edges, narrow out the incisal material thinly towards the shoulder, apply a very thin layer to the tooth body.

D4: Modify dentin D3 with khaki/blue mixture. Only on the edges, draw the incisal material all the way down to the shoulder area, thinly cover tooth body with incisal material.

out after 2 mm (see table) and the desired dentin shape must be completed in a second step.

- Complete shaping with help of the desired **incisal materials** or the effect materials. Due to their transparency and coloring, the **effect materials** are also ideal for individualizing and characterizing the veneers.

Production of fixed inlays, onlays, indirect veneers or fully veneered individual tooth restorations without framework support:

- **Note:** The incisal and effect masses must **not** be applied **directly** to a framework, a bonder surface or an opaque surface. In a first step, a layer of dentin material is therefore applied to the framework after the bonding agent and/or opaquer has been applied. This serves as an elastic or flexible intermediate layer and absorbs forces and stresses. After layering the incisal materials or effect materials, the final polymerization is carried out, see polymerization table.

Curing / Polymerization

- Carry out the final polymerization of a restoration in a suitable light polymerization unit (see polymerization table).
- Depending on the light conditions, the working range of the masses is 1-3 min.
- A suitable light polymerization device with an emission spectrum in the range of at least 310-500 nm must be used. The required physical properties can only be achieved with perfect lamps. It is therefore necessary to regularly check the light intensity according to the manufacturer's specifications.

Note

The exposure time for intermediate polymerization is approx. 60-90 seconds per layer for all colors using a commercially available polymerization device. The final color and finishing is only achieved when the final polymerization / color finishing (color combination table) is carried out after the last layer.

Surface finishing

- After final polymerization, the material may have an increased "yellow value" due to the catalyst. In order to achieve the final color and fix it, a final coating is required depending on the polymerization device used (see table).

Polymerization times for the overpressing technique, see application)

Polymerization-unit	dialog Vario / dialog Vario Flow			dialog Vario Occlusal			
	Paste Opaque	Intermediate polymerization	Final polymerization	Surface finishing	Intermediate polymerization	Final polymerization	Surface finishing
Spektra LED	1 min.	30 sec.	3 min.	none	30 sec.	3 min.	none
Spektra 2000	3 min.	90 sec.	9 min.	7 min.	90 sec.	9 min.	7 min.
HiLite	90 sec.	90 sec.	180 sec.	180 sec.	90 sec.	180 sec.	180 sec.
Spektramat	2 min.	1 min.	5 min.	5 min.	1 min.	10 min.	10 min.
Labolight LV-II / III	1 min.	0,5 min.	9 min.	9 min.	0,5 min.	9 min.	9 min.
Solidlite EX	1 min.	90 sec.	9 min.	7 min.	90 sec.	9 min.	7 min.

- We recommend covering the entire veneering surface with a covering gel before final polymerization. This prevents the formation of a new dispersion layer, guarantees complete polymerization and facilitates finishing.

Application of the composite in the overpress / flask technique (veneering of crowns, bridges, adhesive / bonding bridges, implant-supported restorations and superstructures, telescopic and conical crowns)

The overpress technique facilitates the time-saving and effective veneering of crowns and bridges as well as the esthetic parts of partial dentures and implant-supported dentures with light-curing composite.

Preparation

- After the framework has been created and prepared, the final wax tooth mold is applied to the restoration.
- The flask is used according to the manufacturer's instructions. After boiling out the wax, you should have three molds: A silicone key to control the desired final shape, a clear silicone key in the overpress flask, and a base with exhaust grooves (to catch any excess) for final placement of your design in the flask.
- Sandblasting as well as the application of bonding agent and opaque are carried out in accordance with the manufacturer's instructions.

Layering and polymerization of the materials in the overpress technique

- Apply the crown and bridge paste dentin in the desired shade. Close the flask halfway and place it in a warm place (e. g. evaporator), with a dark cover.
- Close and screw-fix the flask after approx. 10 - 15 min. and polymerize the material in a suitable light curing device, i. e. a device with a sufficiently large interior to accommodate the overpress flask.
- Polymerization time: 5 min.
- Remove the flask from the light curing device. After opening, remove any excess and lift the work from the silicone base. Place the restoration back into the light-curing unit to polymerize the basal surfaces and undercuts.
- Polymerization time: 2 min.
- Grind down the area which is to be completed with incisal material subsequently (see surface finishing).
- After finishing, it is necessary to apply a bonder with subsequent polymerization to restore the dispersion layer (please follow the manufacturer's separate instructions).
- If characterization of the work is desired, stains can now be applied and polymerized according to the manufacturer's instructions.
- **Note:** The surface must not be touched or ground after polymerization. Contamination with oil, dust or dirt must be avoided. The dentin core should be coated directly, otherwise the perfect adaptation of the incisal materials cannot be guaranteed.

- Now apply dialog Vario Occlusal and/or one of the effect masses. Close the flask halfway and place it in a warm place (e. g. evaporator), with a dark cover. Close and screw-fix the flask after approx. 10 - 15 min. and polymerize the material in a suitable light curing device, i. e. a device with a sufficiently large interior to accommodate the overpress flask.

- Polymerization time: 5 min.

- Remove the flask from the light curing device, open it, remove any excess material and grind the surfaces until the desired shape has been created (see Finishing the material).
- Then polish as described in the section Polishing.

Finishing the material

- Suitable instruments for polishing are silicon polishers (wheel, lentic, roll) as well as hard-metal grinders or diamond grinding tools.
- **Note:** As with all composites, a fine sanding dust is formed when dialog Vario and dialog Vario Flow are sanded. We recommend working over an extraction system and wearing protective gloves.

Polishing

- The material is polished with e.g. goat-hair brushes and soft linon buffs mounted on a handpiece, using a polishing paste (e. g. dialog Vario polish). Careful surface treatment and polishing is a prerequisite for an optimum result and largely prevents the formation of deposits (nicotine, caffeine, etc.) and the associated color impairment.
- **Note:** For this working step, we also recommend working over an extraction system and wearing protective gloves.
- After polishing, the construction can be cleaned with water and dried with oil-free compressed air.

Modifications and Repairs

- The composite to be modified is roughened up to 2 mm above the correction or repair margin, coated with bonding fluid and polymerized in the light-curing unit to create a new dispersion layer. Please observe the separate manufacturer's instructions.
- The desired composite material is then applied and finished as described above.

Characterization

The effect materials and suitable dental stains, for example, can be used for characterization. The shades white (increases the opacity and shade in the area of the papilla or near the bone), red (increases the intensity of the pink composite and imitates veins) and blue (enhances the color effect in the shapes and imitates smaller veins) are particularly suitable.

Processing range of the masses

Depending on the light conditions 1 - 3 min.

Troubleshooting

Error	Cause	Remedy
composite does not cure	applied layer of material is too thick insufficient polymerization	maximum layer thickness 2 mm observe polymerization times check lamp; replace lamp if necessary
Layered material chips	opaque not polymerized correctly/ applied layer is too thick faulty framework design	opaquer only applied "wash-like" or in very thin layers incorrect polymerization times / check lamp/replace lamp if necessary
Air bubbles	masses mixed	optimal framework design to ensure that the masses are supported (protection of the mastigratory edges) avoidance of premature contacts do not mix masses, but layer them on top of each other
Mass too dense inside the syringe	syringe was stored at a temperature below 10 °C	store masses at room temperature or let them reach room temperature before
Surface remains greasy	insufficient polymerization	observe polymerization times check unit/ regular maintenance

Material cannot be polished	insufficient polymerization matt surface	observe polymerization times check lamp / replace lamp, if necessary polish sufficiently use suitable polishing paste
Discolorations and increased plaque deposits	insufficient polymerization	check the light output of the polymerization system correct positioning of the objects inside the unit observe the polymerization times of the composites observe the maximum layer strength of the masses seal the surface with a good polish
Air bubbles	insufficient surface coating faulty layering	do not remove the paste directly from the syringe with the instrument, push out and remove the material apply sufficient mass and spread do not mix masses, but layer them on top of each other

Polymerization

- The polymerization depth is 2 mm. For optimum physical properties, a layer thickness of 2 mm should therefore not be exceeded.
- Use a light curing device with an emission spectrum in the range of 310 - 200 nm. The required physical properties can only be achieved with perfect lamps. For this reason, the light intensity must be checked regularly in accordance with the manufacturer's specifications.

Storage

Store the material at 10 - 25 °C to maintain storage stability. Close the syringe tightly immediately after use and protect from direct exposure to light. Turn the spindle back one turn to prevent unwanted leakage of the material.

Shelf-Life

The maximum shelf-life is printed on the label of each packaging. Do not use after the expiration date.

Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (e. g. allergies) or local discomfort, however, cannot be ruled out completely. All serious incidents occurring in connection with the use of the product must be reported to the manufacturer indicated below and to the respective competent authority.

Note on disposal

Residual quantities and packaging material must be disposed of in accordance with local and / or legal regulations.

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Schütz Dental GmbH