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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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#### Sebond Smart

#### Purpose

Sebond Smart is used to bond various restorative materials (e. g. composites, acrylics) to the framework (e. g. metal, zirconium oxide, acrylic) inside the mouth. It is not intended for direct bonding to the tooth structure.

#### Product description

Sebond Smart is a universal bonding agent for a safe and lasting metal-acrylic bond. Sebond Smart is easy to apply without light-curing on all metals and with all conventional composite systems.

#### Users

For application by dental technicians in a dental lab.

#### Composition

Methyl methacrylate, macromeres with phosphonic acid and sulfide groups.

#### Indication

Sebond Smart serves as a bonding agent between acrylics / composites and metals of all kinds.

#### Contraindications

In case of hypersensitivity of the patient to any of the ingredients, this product must not be used or only under strict supervision of the attending physician / dentist.

#### Interaction with other agents

The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

#### Application as bonding agent between acrylics / composites and metals of all kinds

#### Preparatory tasks

After finishing, the surfaces to be veneered must be sand-blasted with aluminium oxide (50 µm, 2 bar). We recommend blowing off the surface activated in this way with oil-free compressed air. Alternatively, the surface can be cleaned with steam or acetone after sand-blasting. In both these cases, sufficient drying of the metal surface must subsequently be ensured.

#### Application of Sebond Smart

Apply Sebond Smart generously with a brush and let it air vor about 1 min. If necessary, apply a second layer.

#### Note

Please provide all of the above information to the treating dentist if you are processing this medical device for custom fabrication. When processing, observe the safety data sheets that exist for this purpose.

#### Troubleshooting/FAQ list

Trouble	Cause	Remedy
bond does not last	applied irregularly or with gaps	pay attention to creating an even film, if necessary apply a second layer
	impurities	keep surfaces clean and grease-free
material evaporated, dried up	bottle was left open	close bottle after each use

#### Storage and handling instructions

Store the material at a temperature between 10° C and 25° C (50° F and 77° F). Close bottle tightly after each use. Protect from direct light incidence.

#### Shelf-life

The maximum shelf-life is printed on the label of each packaging. Do not use after the expiration date.

#### Note on side effects

With proper use of this medical device, unwanted side-effects are extremely rare. All serious incidents occurring in connection with the use of the product must be reported to the manufacturer indicated below and to the respective competent authority.

#### Note on disposal

Residual quantities and packaging material must be disposed of in accordance with local and/or statutory regulations.

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Schütz Dental GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany