



Since 1963

# IMPLA

## No Lock Ti-Base Abutment



### No Lock Ti-Base Abutment

One of the latest innovations in the field of No Lock Ti-Base abutments is the „two-piece“ abutment design. A customized two-piece abutment offers numerous advantages compared to a one-piece abutment.

- Ideal application area for screwed implants
- Works on Hex and Cone Connection each 3.3 mm, 4.2 mm and 5.3 mm
- The new IMPLA No Lock Ti-Base Abutment supports you with an easy insertion of your implant-supported restorations
- Very low insertion height and maximum stability thanks to the metallic connection



### Sebond Implant & Alphalink Implant

This system was specially developed for bonding titanium abutments to individualized zirconium dioxide abutments and supra-constructions. It fixes supra-constructions safely on individual abutments.

- bonding with a high level of safety
- excellent adhesion
- precise, time-saving

Art.-Nr.	No Lock Ti-Base Abutment
637181	3.3 mm Hex connection
637182	4.2 mm Hex connection
637183	5.3 mm Hex connection
637184	3.3 mm Cone connection
637185	4.2 mm Cone connection
637186	5.3 mm Cone connection
640076	Alphalink Implant 8 g
640075	Sebond Implant 1 x 5 ml

Illustrations may vary. We reserve the right to make changes to the product and technical data.



**Do you need any more information?**  
Then arrange a consultation with our experts!



Phone: +49 6003 814-362  
E-Mail: [export@schuetz-dental.de](mailto:export@schuetz-dental.de)